23

pro**3d**ure
medical
DENTAL

Materials are our **DNA!**

Inhalt

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Contents

- 04 Introduction
- 06 Digital workflow
- 08 **printo**dent® NEWS
- 10 **MSI**® technology
- 14 **printo**dent® 3D printing materials
- 48 Devices (3D printer, cleaning, post curing)
- 62 **PRO3**dent®
- 66 **THERMEO**® – the smart mate
- 74 Accessories

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... is the guiding principle of the dynamic, innovative company **pro3dure medical**. The company's philosophy, focused on meeting specific customer requirements, is based on an innate strategic imperative that has been the driving force in researching and developing innovative solutions for the dental world from the very beginning. In doing so, **pro3dure** combines a wealth of experience and concentrated technical know-how with an innovative spirit in order to meet customer-specific requests in a timely manner. **pro3dure** provides customized and sophisticated technology packages for this purpose. The focus is on generative as well as subtractive manufacturing processes and their embedding in digital process chains. We develop customized material concepts for your digital workflow.

Let our team give you an insight into highly innovative material concepts and at the same time an outlook on future technology packages for the production of dental moldings. In line with our IDS motto this year "Speed, Precision, Handling", our focus is on efficient and validated manufacturing processes based on functional materials. Examples include the integration of **printodent®** materials into the **exocad** software packages, the combination with other products (e.g. **Vita Vionic® Bond**) and the expansion of our 3D printer portfolio to include **Prusa's Medical One System**. Another highlight is the result of the collaboration with **Penrhos Bio**, a subsidiary of Unilever. Using the globally unique **MSI® technology** (Multi Species Inhibition), biomimetic polymers for additive and subtractive processes are being demonstrated for the first time that prevent biofilm adhesion and formation. The effect is based on a natural substance formed by the red alga *Delisea* and, unlike bactericides, for example, does not generate bacterial resistance. Accordingly, our materials are not only constantly improving in terms of mechanical properties, but are also getting smarter to achieve better clinical results. Validated processes in combination with material innovations give you the security to be able to work flexibly and efficiently under the rapidly changing regulatory requirements in the future.

Discuss your material requirements for 3D printing and conventional technologies with us. We support you in the development of your product ideas around the digital workflow and we are available as consultants for the implementation of 3D-based process chains. From biometric data to matter: With selected partner companies and a highly qualified team, which has many years of experience and expertise in material and process development, we answer the question for you: What is an efficient way to get biometric data to and from digital back into a real world product? Our quality management system is the important foundation for this. **pro3dure** is certified according to the international standard ISO 13485 and meets additional requirements of participating MDSAP countries. As a result, we guarantee the highest product quality and high-quality customer service worldwide.

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

World's largest selection of 3D printing resins for dental applications by pro3dure ...

From the largest selection of 3D printing resins on the market, start your digital process chain with us. The pro3dure high-tech materials are perfectly matched to a wide range of 3D printers and post-processing units (see page 58-61). This gives you access to the largest variety of dental applications using 3D printing with validated processes.

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The world's largest portfolio of 3D printing resins for dental applications becomes even more comprehensive. Find all our news quickly and follow the hand. **printodent[®]** resins with the worldwide unique **MSI[®] technology** are especially marked by the orange colored triangle.

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
Bacteria have a fascinating and diverse social life. They exhibit coordinated group behavior controlled by **quorum-sensing systems** that detect the density of other bacteria in their environment. The term quorum dates back to the Roman Empire and referred to the lowest number of members required for a vote in the Senate. Similar to the Roman Senate, bacteria vote their behavior as soon as a certain number of other bacteria are in their vicinity. A key example of such group behavior is biofilm formation (e. g. plaque), in which communities of cells attach to a surface and form a polymer shell that protects the bacteria from external influences.

Communication is everything and represents one of the most fundamental mechanisms of intra- and interspecies interactions. Communication can take place in many different ways. One form of human communication is the acoustic transmission of signals to transfer information. In addition to acoustic communication, visual communication also plays a major role in the animal kingdom. In an impressive way, this can be observed in marine-living cuttlefish (octopus). These visually communicating organisms are not only able to change their skin color and surface so that they are perfectly camouflaged. They are also able to move platelets in the chromatophores (specialized skin cells) to reflect specific patterns of polarized radiation. It is now thought that they have a language via these different skin changes, consisting of over 50 vocabularies, and can even form sentences in their specific sequence. Both acoustic and visual communication require that the sender and receiver have highly complex and specialized organs such as ears and eyes.

But what happens when, as in the case of bacteria, the organism cannot possess these organs because it consists of only a single cell? Although bacteria are only single-celled, they can communicate with conspecifics and bacteria of other species. This mechanism was first described in symbiotic bioluminescent bacteria by Nealson et. al. in 1970 and termed **Quorum Sensing** by Fuqua in 1994 and represents a biochemical language. The vocabulary of this language are so-called autoinducers, which are formed by the bacteria themselves and released into the space surrounding. The more bacteria exist in a spatial proximity to each other, the more of these signal molecules are present in this environment. When the concentration of the messenger molecules reaches a certain threshold value (the quorum), the autoinducers flow back into the bacteria, switch on some genes, switch off others and thereby abruptly change the behavior of the organisms. Sometimes the microorganisms form a biofilm that protects them from antibiotics and disinfectants; mouth germs then form harmful dental plaque. In other cases, they produce toxins and decomposing enzymes.

multi-species inhibition

MSI[®] technology

Inspired by Nature | Addicted to Function

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Prof. Richard Hammond, Penrhos Bio Chief Executive Officer:

"We are thrilled to be working with our colleagues **pro3dure medical** to combine the power of **Remora™** technology with their world leading expertise in material development. The result is a completely novel way of inhibiting the formation of harmful multispecies biofilms in their products. **Remora™** is inspired by the natural biology of certain seaweeds that are able to block communication between microorganisms. The result is a surface that remains clean and free of biofilms. Backed by Unilever, we have proven efficacy in a wide range of applications and we are really excited to bring this benefit into dental markets through the innovation of **pro3dure**, an ideal partnership."

Abb. 2: printo[®]dent 3D-Druckharze mit MSI[®] Technologie.
Fig. 2: printo[®]dent 3D printing resins with MSI[®] technology.

REMORA
CLEANER BY NATURE

The formation of biofilms is therefore a significant medical problem (e. g. caries formation). If the communication of microbes could now be interfered with by means of suitable molecules, bacterial infections could also be contained. Unlike antibiotics, the germs would not be killed, but cut off from their information channels. For this reason, the organisms cannot develop resistance to the interfering molecules, as is known from antibiotic resistance. (Fig. 1) In nature, biofilms are not observed on all moist surfaces. There are species of algae that are virtually free of biofilms. For example, the red alga *Delisea* produces a variety of natural compounds (furanones) that appear to block quorum sensing by bacteria. Inspired by nature, **Unilever / Penrhos Bio** has developed a biomimetic technology based on these natural compounds to combat biofilms (called **Remora™**). Research results also document that the bacteriostatic effect works not only on bacteria but also on yeasts, fungi and algae. Thus, a technology is available that can effectively combat „multi-species“ biofilms.

In close cooperation between **Penrhos Bio** and **pro3dure medical**, this multi-species inhibition (**MSI[®] technology**) was realized for the first time using natural substances in 3D printing materials and other dental materials. In doing so, **pro3dure medical** is taking an innovative step towards functional materials in the further development of its dental materials in accordance with the **MSI[®] guiding principle** "Inspired by nature, committed to function". Come and find out more details about **MSI[®] technology** and see numerous new, functional dental materials ranging from milling blanks and denture relining materials to 3D printing resins.

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pro3dure's 3D printing resins of the **printo[®]dent[®] GR-10 series** are biocompatible resins for the fabrication of splints and implant drill guides. In addition to its exceptional elongation for this class of material, **printo[®]dent[®] GR-10 guide** is also characterized by its high transparency and UV stability. This makes it easy to insert sleeves into the guide and gives you optimum control during the operation. The material is now also available in the **printo[®]dent[®] GR-10 guide | MSI** version as a functional 3D printing resin with **MSI[®] technology**, which significantly reduces biofilm and plaque build-up on surfaces. **MSI[®]** is based on the patented mode of action of special natural compounds that disrupt the communication of bacteria (quorum sensing). Inspired by these natural compounds, which are produced by the red alga *Delisea*, a sustainable approach to preventing biofilms has been found that, unlike antibiotics, for example, does not cause bacterial resistance.

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MSI[®] technology in combination with a new polymer matrix has been realized in the **printo[®]dent[®] GR-10.1 guide | MSI** material. In addition to the advantages of reduced biofilm adhesion, the new formulation results in outstanding mechanical properties such as fracture toughness and fracture work. Implant drill guides made of **printo[®]dent[®] GR-10.1 guide | MSI** meet the highest stability requirements and have a fracture toughness and fracture work that are many times higher than those of their market competitors (Fig. 6/7, page 16). This minimizes the risk of splint fracture and increases clinical safety. Due to its novel formulation, it has been possible to dispense with teratogenic components, thus significantly increasing body compatibility. All materials of the **printo[®]dent[®] GR-10 series** can be sterilized in autoclaves using validated processes (Fig. 4).

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pro3dure's printo[®]dent[®] GR-11 tray materials are biocompatible 3D printing resins for the additive manufacturing of customized impression and functional impression trays. Due to their high reactivity, buildstyles with up to 300 µm Z-height can be achieved to be efficient in the lab. With the new development of **pro3dure's printo[®]dent[®] GR-11.1 tray**, the biocompatibility and mechanical stability of the built objects have been taken to another level. The new formulation does not have any reproduction-toxic labelling, thus documenting the special biocompatibility for users and patients. Furthermore, **pro3dure's printo[®]dent[®] GR-11.1 tray** has an impact strength that is approximately 50 % higher than any competitor product (Fig. 9). Accordingly, impression trays can be removed from the patient's mouth with higher safety.


Abb. 9: Schlagzähigkeit (ISO 180), kJ/m² ***
Fig. 9: Impact strength (ISO 180), kJ/m² ***

printo[®]dent[®] GR-11.1 tray
1 kg Flasche / bottle
orange / orange
Art.-Nr. / item-no.: D1001135

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
pro3dure's printo[®] GR-12 cast material is a residue-free burn-out 3D printing resin that can be used in combination with a variety of investment materials for a wide range of applications. Thanks to the optimized setting of the curing depth, even the finest details can be printed precisely with the material. The high stability of the material in the green compact mode even allows a direct FIT check. Components made of **printo[®] GR-12 cast** can be cleaned quickly and efficiently thanks to its low viscosity. If, despite their high stability, **printo[®] GR-12 cast** objects need to be repaired or augmented, we recommend the associated UV-curable repair gel **LCR-1** from **pro3dure** (see page 74).

printo[®] GR-12 cast
1 kg Flasche / bottle
rot / red
Art.-Nr. / item-no.: D1000737

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¹ Holding time/casting process according to manufacturer's instructions

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- rückstandsfrei ausbrennbar
- hohe Präzision, feinste Details
- hohe Stabilität
- schnelle, effiziente Reinigung
- reparierbar mit **LCR-1**
- residue-free burn-out
- high precision, finest details
- high stability
- fast, efficient cleaning
- repairable with **LCR-1**

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high dimensional stability printo[®] dent GR-13 model

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- hohe Dimensionsstabilität
- optimierte Optik und Haptik
- Implementierung in den digitalen Workflow
- high dimensional stability
- optimized look and feel
- implementable in the digital workflow

printo[®] GR-13 model
1 kg Flasche / bottle
beige / beige
Art.-Nr. / item-no.: D1000352

hohe Dimensions-
stabilität
look
and feel
optimized

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“3D printing with **printo[®]** in your own medical office is the perfect complement to our intraoral scanner.”

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pro3dure's printo[®] GR-13 model is a 3D printing resin for the generative fabrication of dental models. The material is available as **printo[®] GR-13 model** for restorative dental models and as **printo[®] GR-13.1** / **printo[®] GR-13.2 model** for orthodontic applications. All resins are characterized by high dimensional stability and the optimized look and feel for dental applications. Due to colorant and opacity the important features like margin lines are becoming perfectly visual. Moreover, the implementation of 3D printed models into your digital workflow is supported. Dental models from the **GR-13 series** can be used effectively in thermoforming processes due to their thermal stability.

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The new **printo[®]dent GR-13.3 model release** was specially developed for orthodontic applications and eliminates the need for additional release agents when using 2-component powder/liquid resins with 3D-printed dental models. Through **pro3dure**'s new "release" technology, it is possible for the first time to provide a material that has an integrated separation effect. This is not due to the addition of fillers such as waxes, which can segregate from the resin matrix and lead to unstable builds or defective components. Furthermore, due to its low viscosity, components made of **printo[®]dent GR-13.3 model release** are easy to clean, thus enabling an efficient finishing process. It is suitable for many othodontic powder/liquid systems[°].

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

multi-species
inhibition

GR-14.2 denture HI | MSI

a) Bauplattform mit Arbeitsmodell (inklusive separater Gingiva-Masken) aus printo[®] GR-13.3 model release. / Building platform with working model (including separate gingival masks) made of printo[®] GR-13.3 model release.

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Combination of different
model materials
in the digital workflow

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- easy handling
- faster working
- fewer work steps
- no separator necessary

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The printo[®] GR-14 denture materials are 3D printing resins for the generative production of denture bases. The printo[®] GR-14.2 denture HI meets all dental requirements with e.g. a water absorption ($< 25 \mu\text{g}/\text{mm}^3$) and water solubility ($< 1,5 \mu\text{g}/\text{mm}^3$) according to ISO 20795-1:2013. Furthermore, the low shrinkage (approx. $< 1 \%$) compared to traditionally used PMMA-based materials guarantees a high accuracy of fit. printo[®] GR-14.2 denture HI is a 3D printing resin with „high-impact“ properties characterized by low viscosity. This allows construction processes to be carried out efficiently and components to be cleaned easily. As a world first, printo[®] GR-14.2 denture HI | MSI, a biomimetic 3D printing resin for the fabrication of dentures, will be presented for the first time at IDS 2023. This resin offers all the advantages of printo[®] GR-14.2 denture HI, which has been tried and tested thousands of times, and also prevents biofilms from adhering to the surface of the denture in a natural way. This functional material thus copies a natural process, such as that used by the red alga Delisea. The communication of bacteria, the so-called „quorum sensing“, is disturbed by a natural substance anchored in printo[®] GR-14.2 denture HI | MSI, so that the adhesion of bacteria is significantly reduced. Moreover, this effect is not due to antibacterial additives, which may lead to bacterial resistance. As a result, improved oral hygiene and less irritation of the soft tissue are achieved.

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- high accuracy of fit
- „high-impact“ material
- validated processes
- MSI[®] technology for the prevention of biofilms

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Fig. 14 (a-e): Use of different model materials (printo[®] GR-13 model and 13.3 model release) in the digital workflow.

multi-species inhibition printo[®] dent[®] GR-14.3 denture HI | MSI

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- high accuracy of fit
- „high-impact“ material
- validated processes
- TPO-free
- **MSI[®] technology** to prevent biofilms

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“Biofilm inhibition with the Remora[™] brand technology results in cleaner and healthier dental products and appliances which assist in maintaining oral health and stability.”

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As a further member of the **printo[®] dent[®] GR-14 denture** family, the **printo[®] dent[®] GR-14.3 denture HI | MSI** is also equipped with **MSI[®] technology** to prevent biofilm formation. In combination with outstanding chemical-physical properties, dentures made of this material have a thermal effect optimized for body temperature, which generates an extremely high level of wearer comfort. All materials of the **printo[®] dent[®] GR-14 denture series** can of course be combined with **printo[®] dent[®] GR-17.1 temporary It** to produce cost-effective and highly esthetic full dentures. In addition, validated procedures are available for both materials in combination with denture teeth from the **VITA VIONIC VIGO**, **CANDULOR TCR* / NFC**** and **Ruthinium ACRYSMART** systems. Complete integration in **exocad's** CAD software rounds off the digital workflow.

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Fig. 15-17: Validated process for the production of full dentures based on pro3dure's printodent[®] GR-14.2 denture HI in combination with VITA VIONIC VIGO[®] – digital tooth solution.

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pro3dure's printodent[®] GR-15 gingiva is a 3D printing resin for the production of flexible gingiva masks. Its natural-looking gingival color combined with a gingival-like flexibility allow optimal dental work with a view to "red and white aesthetics". The printodent[®] GR-15 gingiva is the perfect complement to pro3dure's printodent[®] GR-13 model to fabricate implant models based on biometric data.

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"To ensure a harmonious transition from the crown to the gum, printodent[®] GR-15 gingiva is excellent."

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printc[®] dent

The largest portfolio of 3D printing
resins for dental applications by ...

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"We print a radiopaque scanning appliance from **printodent[®] GR-16 Xray** with each of our end of the line digital dentures and All-On-X cases. This allows us to distinctly differentiate our lab with our digital workflow and deliver a better prostheses to our dentists."

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pro3dure's printodent[®] GR-16 Xray is the first radiopaque 3D printing resin for the generative fabrication of scanning templates. On the basis of a stabilization technology developed by **pro3dure** for filled materials and a suitable selection of radiopaque additives, the construction of scanning templates can now be integrated into your "digital workflow" efficiently.

Abb./Fig. 19:
Röntgensichtbarkeit **printodent[®] GR-16 Xray** im Vergleich zu einer Aluminiumtreppe (ISO 4049)./
Radiopacity of **printodent[®] GR-16 Xray** in comparison to aluminum stairs (ISO 4049).

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printo[®]dent

Langzeit-
provisorien
longterm
temporaries

By courtesy of:
Prof. Alessandro Pozzi
DDS Ph. D.
Rome, Italy

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
- biocompatible materials
- very good mechanical properties
- double indication
- 7 colors - 2 sizes
- low viscosity
- easy to clean
- radiopaque

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pro3dure's GR-17 series resins are biocompatible materials for the additive manufacturing of temporary crowns and bridges. The **printodent[®] GR-17 temporary** is indicated for temporary restorations in the anterior region. **printodent[®] GR-17.1 temporary It** has exceptional mechanical properties compared to conventional temporary materials with a Flexural modulus > 5000 MPa (according to ISO 178**) and a Flexural strength of 112 MPa (ISO 4049). Moreover **printodent[®] GR-17.1 temporary It** exhibits excellent tissue compatibility and can therefore be used ideally for longterm temporaries¹. As the first material in its class worldwide, **printodent[®] GR-17.1 temporary It** received FDA approval for two indications: temporaries and denture teeth. This expands the possibilities in your laboratory and makes the use of this material extremely efficient and economical. Both materials are available in tooth shades A1, A2, A3, A3.5, B1, B2 and bleach.

¹ Lankes V, Reymus M, Liebermann A, Stawarczyk B. Bond strength between temporary 3D printable resin and conventional resin composite: influence of cleaning methods and air-abrasion parameters. Clin Oral Investig. 2023 Jan;27(1):31-43. doi: 10.1007/s00784-022-04800-7

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Radiopacity of printo[®] GR-17.1
temporary It in comparison to
aluminum stairs (ISO 4049).

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For the temporary restoration after
implantation, I found the optimal
solution for my digital workflow
with printo[®] GR-17.1 temporary It.

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Scientific studies¹ prove the compatibility with bredent
visio.link as well as the optimal adhesive bond with
printo[®] GR-17.1 temporary It.

¹Lankes V, Reymus M, Liebermann A, Stawarczyk B. Bond strength between temporary 3D prin-
table resin and conventional resin composite: influence of cleaning methods and air-abrasion
parameters. Clin Oral Investig. 2023 Jan;27(1):31–43. doi: 10.1007/s00784-022-04800-7

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pro3dure's printodent[®] GR-18.2 IB is a biocompatible 3D printing resin for the additive manufacturing of orthodontically used „Indirect Bonding Trays“. The material is flexible, transparent and has a high tensile strength. This enables the orthodontist to place all brackets quickly and precisely in a single operation. The further development of the proven printodent[®] GR-18.1 IB focused in particular on optimizing application technology. For example, printodent[®] GR-18.2 IB has a lower viscosity and is therefore easier to clean. Furthermore, the inhibition layer has been reduced during post-curing with low-performance post-curing systems. As a result, you avoid sticky component surfaces and achieve higher clinical safety.

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- biocompatible
- flexible and transparent
- high tear resistance
- precise placement
- higher clinical safety

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pro3dure's printodent[®] GR-19.1 OA | MSI is the first biomimetic 3D printing resin for the generative production of hard-elastic orthodontic splints. Thanks to the unique MSI[®] technology, splints made of printodent[®] GR-19.1 OA | MSI exhibit reduced biofilm adhesion. The effect, which is copied from nature, is based on natural substances that sustainably inhibit the communication of bacteria. Unlike antibiotics, for example, no resistance is formed in this way. The material also has a body temperature-induced thermo effect and thus ensures maximum wearing comfort in the mouth.

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- reduced biofilm adhesion due to MSI[®] technology
- hard flexible
- highest wearing comfort due to thermo effect

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pro3dure's printo[®] dent GR-20 MJF is a biocompatible and hemocompatible 3D printing resin for the generative fabrication of oral and maxillofacial surgery devices. printo[®] dent GR-20 MJF parts can be sterilized using plasma, gamma rays, ethylene oxide or an autoclave. Moreover, the material is resilient against disinfectants. As the first of its kind it offers a broad range of applications in every operation room.

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printo[®]dent[®] GR-21.1 Try-In
1 kg Flasche / bottle

A1 Art.-Nr. / item-no.: **D1001651**
A2 Art.-Nr. / item-no.: **D1001652**
A3 Art.-Nr. / item-no.: **D1001653**

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printo[®]dent[®] GR-21.1 Try-In is a class 1 biocompatible medical device and is suitable as base material for functional try-in of digitally fabricated denture bases. **printo[®]dent[®] GR-21.1 Try-In** can be used to generatively fabricate dental moldings such as denture bases with individual tooth set-up in order to perform a functional try-in to check bite registration and occlusion. The material is available in tooth shades A1, A2 and A3. The newly developed **printo[®]dent[®] GR-21.1 Try-In** is characterized by significantly improved chemical and physical values compared to its market competitors. For example, the impact strength (according to ISO 180/U) has been increased by 75 % compared to its predecessor. The use of a novel initiator system has further increased biocompatibility and thus clinical safety.

a) **printo[®]dent[®] GR-21.1 Try-In** gedrucktes Gerüst für die Anprobe einer verschraubbaren Implantatbrücke./ **printo[®]dent[®] GR-21.1 Try-In** printed framework for the try-in of a screw-retained implant bridge.

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Fig. 22 (a-d): Try-in on an implant bridge with **printo[®]dent[®] GR-21.1 Try-In**

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"It has long been my wish to be able to produce flexible splints without impressions. With **printo[®] dent GR-22 flex**, we can now print our flexible splints in high quality. A significant step forward on my digital path."

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printo[®] dent GR-22 flex – the innovation for 3D printing.

The number of treatments of occlusal disorders by means of splint therapy is constantly increasing and thus the importance of the "dental splint" for a wide variety of indications such as bruxism, bite elevations or mouth protection. The demands placed on splint materials by patients, practitioners and dental technicians are manifold and range from durability, esthetics, wearing comfort and biocompatibility to processability and usability in digital processes. Based on biometric data, 3D printing processes can nowadays be used to produce complexly designed therapeutic splints efficiently. The latest member of the printo[®] dent family was developed with these aspects in mind – the innovative, flexible splint material printo[®] dent GR-22 flex with thermomemory effect, optimized for additive manufacturing. After all, who doesn't know it, the challenge of a "crooked impression" during conventional impression taking or imperfections during intraoral scanning that lead to "data holes" on the surface of a generated data set? A conventional, hard splint then wobbles on the plaster model and in the patient's mouth. The splint is subjected to greater stress at precisely these points. The risk of chipping or breaking is significantly increased. In addition, the wearing comfort for the patient is reduced. The innovative printo[®] dent GR-22 flex compensates for such inaccuracies thanks to its elastic properties combined with a thermomemory effect. printo[®] dent GR-22 flex becomes flexible at body temperature and thus adapts optimally to the tooth surface. This increases wearing comfort and minimizes the risk of fracture. printo[®] dent GR-22 flex thus supports the satisfaction of dental technicians, practitioners and patients through higher safety, exceptional wearing comfort and fewer "recalls".

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printodent® GR-22 flex is the flexible 3D printing resin for splints.

printodent® GR-22 flex gives you and your customers the flexibility you need! The development of the printodent® GR-22 flex was based on extensive clinical and material data from conventional, flexible splint materials. As a result, the printodent® GR-22 flex achieves increased flexibility in combination with a more pronounced thermomemory effect compared to its market competitors (Fig. 23), which is immediately noticeable when the splint is inserted in the mouth due to a higher pressure drop at the tooth. The wearing comfort for the patient is thus significantly increased. Would you also like to remain flexible in the selection and design of your "digital workflow"? A variety of validated processes based on a wide range of CAD softwares, 3D printers and post-treatment units are available for the printodent® GR-22 flex.

inside **exocad**

be outstanding

printodent® GR-22 flex – a resin with outstanding properties.

printodent® GR 22 flex meets the highest medical technology standards. Dental splints made of printodent® GR 22 flex have outstanding mechanical properties compared to their market competitors. For example, the standard specifications for impact strength (ASTM_D256) and elongation at break (ISO 37) are not only achieved, but exceeded by more than 3.5 and 2.5 times respectively (Fig. 24 and Fig. 25). The new development printodent® GR 22 flex is based on an innovative polymer matrix which ensures that the solubility of polymer components (Fig. 26) is also significantly lower than that of the comparative products. This ensures outstanding biocompatibility in this class of materials.

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printodent® GR-22 flex is durable and biocompatible.

Characterized by its exceptional wearing comfort, printodent® GR-22 flex has been developed for durability and reliability. The stability of the material was proven in an extensive chewing simulation study conducted by LMU Munich¹. Although printodent® GR-22 flex has a much higher flexibility than its market competitor, no significant differences in terms of material loss (volumetric and vertical) could be found between them (Fig. 27 and Fig. 28). printodent® GR-22 flex thus increases patient comfort without reducing clinical benefit! Furthermore, tested reworking and cleaning recommendations are available to sustainably ensure the special properties of a printodent® GR-22 flex splint over its entire service life.

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Fig. 27: printodent® GR-22 flex (orange) and KeySplint Soft (blue) volumetric degradation after 120,000 cycles chewing simulation, mm³.¹

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Fig. 28: printodent® GR-22 flex (orange) and KeySplint Soft (blue) vertical degradation after 120,000 cycles chewing simulation, mm.¹

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printodent® GR-22 flex is the innovative resin for 3D printing of splints with unique impact strength combined with high elongation at break. Its special viscoplastic properties in combination with a thermomemory effect combine highest comfort with highest biocompatibility.

- high elongation at break
- unique impact strength
- highest wearing comfort (usability studies pro3dure)
- lowest solubility/high biocompatibility
- thermomemory effect
- highest impact strength in its class
- "snap on" splints available

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Compact, user-friendly and precise are the attributes of the **Prusa Medical Systems**, which consists of the **Prusa Medical One** 3D printer and the **CW One** washing and post-curing unit. The Prusa Medical One has a high-resolution UV-LED panel and belongs to the group of image projection systems based on LED-LCD technology. A variety of special solutions, such as build platforms specially adapted to dental requirements and software features to compensate for curing layers, guarantee precise and reproducible results in everyday work. In combination with the all-in-one solution CW One, which allows you to wash and post-cure your components, a complete process chain for the production of dental molded parts can be realized in a minimum of space at an exceptional price level. With this system, you have access to validated processes with **pro3dure's printodent** materials and thus to the world's largest portfolio of dental applications for 3D printing. The indications range from precise dental models with **printodent® GR-13 model**, splints and drill templates with **printodent® GR-10.1guide | MSI**, denture bases with **printodent® GR 14.2 HI | MSI**, denture teeth, crowns and bridges for temporary long-term care with **printodent® GR-17.1 temporary It**. Prusa Medical One convinces with precise printing results and easy handling.

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"With **pro3dure** and Prusa, I have all validated parameters directly from the manufacturer. So my workflow is precise and accurate."

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"We are very happy that we have teamed up with **pro3dure**! They provide top-notch biocompatible resins and are excellent business partners."

Josef Prusa
(CEO and Founder Prusa Research)

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Fig. 29: Precision printing platforms are an important part of our Medical One 3D printer. They are manufactured with maximum precision so that they can be freely exchanged without the need to calibrate the printer. This results in a significant reduction of the downtime between prints and even higher efficiency of the entire printing process.

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Fig. 30: The **Prusa Slicer** software package features a large number of special solutions that are specifically tailored to the requirements of the dental sector.

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The combination of the **pro3dure** materials and the proven **ASIGA MAX™** form the backbone of your digital process chain. The ASIGA MAX™ is the world's most advanced digital dental 3D printer. It features maximum productivity while requiring a minimal physical footprint. With 62 µm, HD print precision and a 385 nm UV irradiation source, the complete portfolio of **pro3dure's** materials can be used. In order to get access to the largest dental application spectrum currently on the market, **printodont®** by **pro3dure**, the ASIGA MAX™ is your "All-round" 3D printer.

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The **Asiga PRO family** combines proven industry leading precision, reliability and large build envelopes which deliver high-speed print capability for the professional environment. You have the option to choose between several lateral resolutions and build envelopes to fit your needs. The new PRO family can be integrated perfectly into your own production process delivering the highest quality and efficiency.

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
Processor: 2.4 GHz multi-core
Memory: 4 GB
Graphic: 256 MB, OpenGL 2.0 or later support
Mouse: 3-button mouse with scroll wheel
Network: Required for printing
Storage: 1 GB free

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Enhance the performance of your **Asiga Max™** system in your dental environment by using the **Asiga Composer software package**. Asiga Composer acts as a link between your CAD-Software and your Asiga printer in order to realize a smooth and integrated digital workflow. The Composer is an intuitive, automatic STL & SLC support and build preparation package. Based on the novel optimal support algorithm, parts are supported automatically only where absolutely required by the part geometry. This makes job setup and post-processing a breeze. With Composer you can include any STL, SLC and Stomp files in the same build. This saves users from having to convert files from one type to another. Projects can be saved as a lightweight "Asiga Composer Build" (.comblb) or complete export into a ZIP file. Composer's automatic part placement tool handles efficient placement of both STL & SLC files based on part contours. This maximizes your available build area. You can manage your workflow by preparing more than one print at the same time and switch between platforms which are conveniently tabbed at the top of your screen. Composer is a cross platform desktop software supporting Windows, Macintosh and Linux. It is included free with every system with an unlimited seats license for use throughout your office.

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The **Miicraft Prime and Profession series printers** offer a wide range of build platform sizes, wavelengths, and xy resolutions to efficiently meet a wide variety of dental applications. In addition, the proven Miicraft beamer technology gives you the confidence to cope with changing demands in the daily laboratory routine. Of course, validated manufacturing processes for these 3D printers are available for all **printodent®** materials.

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The new "cleaning device" **CLD-1** represents another building block for the optimized and controlled manufacturing process of medical devices by generative fabrication. Unlike conventional ultrasonic baths the **CLD-1** generates a vortex of the cleaning agent. As a result, microcracks in the component surfaces are significantly reduced, thus further increasing the part quality in terms of surface finish, appearance and mechanics of your dental objects. A selection of different cleaning programs assists you in perfectly matching the cleaning program to the selected material class and choice of cleaning agent. The cleaning container can be easily removed and changed in order to simplify the manufacturing processes in your lab. Changing the cleaning liquid can be done within seconds in a convenient way.

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The **RLR-1** is a roller/tilting stirring device for homogenization of 3D printing materials with pigment or filler content. Highly filled 3D printing materials like the **printodent® GR-17.1 temporary It** require a thorough mixing before use, handshaking is often insufficient. This can result in significant colour deviations and even build failures may occur. The RLR-1 Mixer is constructed for constant operation in order to keep your **pro3dure** 3D materials ready for use at any time with the same consistency.

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

The high-tech post-curing unit **CD-2** is tuned to the latest generation of 3D printers. It features top-level productivity, the most current LED irradiation technology, sophisticated sensors and a wide range of smart features. Through the touch panel you can fully control and monitor all relevant irradiation parameters. Sensors measure the protective gas flow and inform the user about the status of the light source. This avoids inhibition layers on generated parts. As a result, the biocompatibility and the optimal mechanical properties of your 3D printed parts are guaranteed. In a profound and unprecedented form the CD-2 combines compactness and efficiency on the highest level of technology, along with an excellent price/performance ratio.

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printo[®]
dent

		printo[®]dent																			
		GR-10 guideGR-10 guide MSIGR-10.1 guide MSIGR-11.1 trayGR-12 castGR-13 modelGR-13.1 modelGR-13.2 modelGR-13.3 model releaseGR-14.2 denture HI MSIGR-14.3 denture HI MSIGR-15 gingivaGR-16 XrayGR-17 temporaryGR-17.1 temporary ItGR-18.2 IBGR-19.1 OA MSIGR-20 MJFGR-21.1 Try-InGR-22 flex																			
ASIGA MAX		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	
ASIGA MAX UV		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
ASIGA PRO 2		✓	○	○	○	✓	✓	✓	✓	✓	○	○	✓	✓	✓	✓	○	○	✓	○	
ASIGA PRO 4K		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
SprintRay Pro 95		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	
MiiCraft Ultra		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
MiiCraft Prime		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
MiiCraft Profession 4K		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
MiiCraft Alpha		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Organical X1S		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Organical X1N		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
phrozen Shuffle/Sonic		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	

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printo[®]
dent

		GR-10 guide GR-10 guide MSI GR-10.1 guide MSI GR-11.1 tray GR-12 cast GR-13 model GR-13.1 model GR-13.2 model GR-13.3 model release GR-14.2 denture HI GR-14.2 denture HI MSI GR-14.3 denture HI MSI GR-15 gingiva GR-16 Xray GR-17 temporary GR-17.1 temporary It GR-18.2 IB GR-19.1 OA MSI GR-20 MJF GR-21.1 Try-In GR-22 flex																			
SISMA EVERES UNO		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Prusa MEDICAL ONE		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Anycubic Photon Mono X		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Microlay Versus		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Microlay Eve pro		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
SHINING 3D AccuFab-L4D		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
SHINING 3D AccuFab-D1s		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	
Flashforge Focus 6K		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Structo ST-01		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Uniz SLASH PLUS		○	✓	✓	○	○	○	○	○	○	✓	✓	○	○	✓	○	✓	✓	○	○	
UNIZ NBEE		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
NOBIL METAL Sinergia Print		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	

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printo
dent

		GR-10 guide	GR-10 guide MSI	GR-10.1 guide MSI	GR-11.1 tray	GR-12 cast	GR-13 model	GR-13.1 model	GR-13.2 model	GR-13.3 model release	GR-14.2 denture HI	GR-14.2 denture HI MSI	GR-14.3 denture HI MSI	GR-15 gingiva	GR-16 Xray	GR-17 temporary	GR-17.1 temporary It	GR-18.2 IB	GR-19.1 OA MSI	GR-20 MJF	GR-21.1 Try-In	GR-22 flex
pro3dure CD-1		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
pro3dure CD-2		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
NK-OPTIK Otoflash G171		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
ASIGA Flash		○	○	○	○	✓	✓	✓	✓	○	○	○	○	○	○	○	○	○	○	○	○	○
Kulzer cara Print LEDCure		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	○	✓	○	✓	✓	✓	✓
SHINING 3D FabCure		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	○	✓	✓	✓	✓
SprintRay ProCure		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
SprintRay ProCure 2		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
SISMA EVERES CURO		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Prusa MEDICAL CW ONE		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
WICKED ENGINEERING CUREbox		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
rapidshape RS Cure		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Nobil Metal Dentafarm Photopol		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

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printo
dent

		GR-10 guide	GR-10 guide MSI	GR-10.1 guide MSI	GR-11.1 tray	GR-12 cast	GR-13 model	GR-13.1 model	GR-13.2 model	GR-13.3 model release	GR-14.2 denture HI	GR-14.2 denture HI MSI	GR-14.3 denture HI MSI	GR-15 gingiva	GR-16 Xray	GR-17 temporary	GR-17.1 temporary It	GR-18.2 IB	GR-19.1 OA MSI	GR-20 MJF	GR-21.1 Try-In	GR-22 flex
pro3dure CLD-1		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
SprintRay Pro Wash/Dry		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Kulzer cara Print Clean pro		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Prusa MEDICAL CW ONE		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
rapidshape RS Wash		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

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To complement our 3D printing resin portfolio **printodent[®]** and the high-tech **THERMEO[®]** resins, we bundle resins for a wide range of prosthetic applications under the umbrella brand **PRO3dent[®]**. For the digital and analog workflow, you will find technically outstanding products with the latest technology under the **PRO3dent[®]** brand. Optionally equipped with **MSI[®] technology**, **PRO3dent[®]** products offer a unique additional feature: multi-species inhibition significantly reduces biofilm accumulation and prevents plaque build-up.

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PRO3dent[®] splint | MSI is ideal for CAM-fabricated hard bite splints. Multi-species inhibition significantly reduces biofilm accumulation and prevents plaque build-up. Available in ocean-blue, in thicknesses of 16 mm and 20 mm (single packaging Ronde 98.5 mm with step).

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- hohe Biegefestigkeit
- keine Verfärbungen
- höchste Biokompatibilität
- Multi-Species-Inhibierung (**MSI[®]**)
- signifikante Reduktion von Biofilmen
- geringe Plaqueanlagerung
- high flexural strength
- no discoloration
- highest biocompatibility
- multi-species inhibition (**MSI[®]**)
- significant reduction of biofilms
- low plaque accumulation

PRO3dent [®] splint MSI	Ronde ozean-blau / blank ocean-blue	16 mm	Art.-Nr. / item-no.: D1025002
PRO3dent [®] splint MSI	Ronde ozean-blau / blank ocean-blue	20 mm	Art.-Nr. / item-no.: D1025003

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PRO3dent[®] base | MSI is used for the CAM production of individual denture bases in the digital workflow. Available in three colors and two thicknesses, the round blank with a diameter of 98.5 mm offers outstanding properties. The high flexural strength as well as the excellent bending modulus ensure optimum stability. Multi-species inhibition significantly reduces biofilm accumulation and prevents plaque build-up. Available in light pink, orange-pink and dark pink, in thicknesses of 25 mm and 30 mm (single packaging round 98.5 mm with step).

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- high flexural strength
- no discoloration
- highest biocompatibility
- multi-species inhibition (MSI[®])
- significant reduction of biofilms
- low plaque accumulation

PRO3dent [®] base MSI [®]	Ronde hell-pink / blank light-pink	25 mm	Art.-Nr. / item-no.: D1024016
PRO3dent [®] base MSI [®]	Ronde hell-pink / blank light-pink	30 mm	Art.-Nr. / item-no.: D1024017
PRO3dent [®] base MSI [®]	Ronde orange-pink / blank orange-pink	25 mm	Art.-Nr. / item-no.: D1024026
PRO3dent [®] base MSI [®]	Ronde orange-pink / blank orange-pink	30 mm	Art.-Nr. / item-no.: D1024027
PRO3dent [®] base MSI [®]	Ronde dunkel-pink / blank deep-pink	25 mm	Art.-Nr. / item-no.: D1024006
PRO3dent [®] base MSI [®]	Ronde dunkel-pink / blank deep-pink	30 mm	Art.-Nr. / item-no.: D1024007

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... is a permanently hard-flexible relining material on a novel polymer base for total, partial and implant-anchored dentures. **PRO3dent[®] reliner** has convincing advantages such as an outstanding adhesive bond with all commercially available denture materials and a tear and tear propagation resistance that is unsurpassed in its class. In addition, its silk-matt surface ensures an outstanding suction effect of the prosthesis in combination with highest wearing comfort. In addition to its exceptional mechanical properties, **PRO3dent[®] reliner** is also the world's first biomimetic relining material that naturally avoids biofilms in the oral environment through **pro3dure**s unique **MSI[®] technology**.

PRO3dent[®] reliner | MSI[®]
75 g Dose / can
hell-pink / light-pink
Art.-Nr. / item-no.: D1001801

PRO3dent[®] reliner | MSI[®]
75 g Dose / can
orange-pink / orange-pink
Art.-Nr. / item-no.: D1001802

PRO3dent[®] reliner | MSI[®]
75 g Dose / can
dunkel-pink / deep-pink
Art.-Nr. / item-no.: D1001800

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THERMEO® is flexible.

Who doesn't know it, "crooked impressions" during conventional impression taking or imperfections during intraoral scanning that lead to "data holes" on the surface of the generated data set? A conventional hard splint then wobbles on the plaster model and thus in the patient's mouth. The splint is subjected to greater stress at this point, the risk of fracture is significantly increased and the patient complains about a lack of wearing comfort. The unique THERMEO® splint material compensates for such inaccuracies with its patented thermomemory effect. The THERMEO® material becomes flexible at body temperature and adapts optimally to the tooth surface. This increases the wearing comfort of the splint and minimizes the risk of breakage. THERMEO® thus supports the satisfaction of the dental technician, dentist and patient through higher safety, exceptional wearing comfort and fewer "recalls".

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"We are proud to be able to equip internationally successful world-class athletes in the field of sports dentistry with THERMEO® splints with thermomemory effect."

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THERMEO® is smart.

Smart materials/responsive materials are materials that react independently to changing environmental conditions such as temperature. This reaction is reversible. The high-tech plastics of the THERMEO® family belong to this group of smart materials and are characterized by a unique patented thermomemory effect. Compared to hard dental splints, inaccuracies during impression taking can thus be compensated with a THERMEO® splint. THERMEO® products can also be used to produce functional therapeutic bite splints. THERMEO® splints are characterized by high wearing comfort and are self-adjusting due to the thermomemory effect. THERMEO®s special properties combine the advantages of hard as well as soft splint materials.

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Fig. 31: THERMEO® (orange) and competitor A (red) vertical degradation after 120,000 cycles chewing simulation, mm.*

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Fig. 32: THERMEO® (orange) and competitor A (red) volumetric degradation after 120,000 cycles chewing simulation, mm³.*

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THERMEO® is durable.

The materials of the THERMEO® family, characterized by their exceptional wearing comfort, have been developed for durability. The stability of the material was documented in an extensive chewing simulation study conducted by LMU Munich*. Despite the higher flexibility and higher resilience of THERMEO® compared to its market competitors, comparable abrasion values were determined (Fig. 31 and 32). Thus, patient comfort is increased without reducing the clinical benefit. Furthermore, amine compounds and MMA were avoided in the development of THERMEO® in order to reduce the yellowing of the material and its potential for irritation and sensitization. Furthermore, tested cleaning recommendations for THERMEO® are available to the patient in order to sustainably ensure the special properties over the entire service life of the THERMEO® splint (see page 73).

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Fig. 34: THERMEO® time depending stress relaxation at 73.4 °F (23 °C).*

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Fig. 35: THERMEO® (orange) and competitor A (red) decrease in stress at 98.6 °F (37 °C) in comparison to 73.4 °F (23 °C).*

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Fig. 36: THERMEO® (orange) and competitor A (red) elastic recovery after a 90° deformation depending on the temperature, 73.4 °F (23 °C) and 98.6 °F (37 °C).*

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THERMEO® is comfortable.

THERMEO® meets the highest medical technology standards. For example, THERMEO® does not use phthalate-based plasticizers, but instead employs a patented flexibilization technology that ensures unique biocompatibility in this class of material. Compared to other suppliers, this results in a significantly greater pressure drop on the tooth after insertion of the splint (Fig. 34 and 35), which in turn increases wearing comfort for the patient. Furthermore, due to the technology used, the thermomemory effect and thus the resilience of the THERMEO® material is greater than that of competitors (Fig. 36). The THERMEO® splint returns to its original position faster and more precisely due to body temperature. THERMEO® can be processed conventionally using the casting or tamping-pressing technique or integrated into the digital workflow in the form of THERMEO® blanks. The blanks are available in thicknesses of 16, 20 and 25 mm. This allows tooth-supported and also soft tissue-supported splints to be manufactured efficiently. Both the powder/liquid system and the milling rounds are additionally available as THERMEO® SO in shades A1, A2 and A3 for the fabrication of "snap-on splints" (Fig. 33).

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The thousands of times proven THERMEO® round is now also available as an MSI® (Multi Species Inhibition) version. This makes it the world's first biomimetic milling disc in its class, which naturally prevents the adhesion of biofilms in the oral environment using pro3dure's unique MSI® technology. It copies a natural process, such as that used by the red alga Delisea. A natural substance (lactam) anchored in the dental material interferes with the communication of bacteria, the so-called „quorum sensing“. This significantly reduces the formation of biofilms or their adhesion. This effect is not due to antibacterial additives, which may lead to bacterial resistance. In the course of this, improved oral hygiene is achieved and thus less irritation of the soft tissue in accordance with the motto „Inspired by nature, addicted to function“.

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- exceptional wearing comfort
- unique flexibility by thermomemory effect
- self adjusting
- no yellowing
- extreme durability
- long-lasting
- highest biocompatibility
- MSI® technology for the prevention of biofilms

multi-species inhibition inside exocad

THERMEO®	Starterkit / Starterkit		Art.-Nr. / item-no.: D1010118
THERMEO®	Flüssigkeit / liquid	130 ml	Art.-Nr. / item-no.: D1010122
THERMEO®	Flüssigkeit / liquid	255 ml	Art.-Nr. / item-no.: D1010123
THERMEO®	Flüssigkeit / liquid	600 ml	Art.-Nr. / item-no.: D1010124
THERMEO®	Pulver / powder	180 g	Art.-Nr. / item-no.: D1010119
THERMEO®	Pulver / powder	340 g	Art.-Nr. / item-no.: D1010120
THERMEO®	Pulver / powder	800 g	Art.-Nr. / item-no.: D1010121
THERMEO®	Ronde / blank	16 mm	Art.-Nr. / item-no.: D1022002
THERMEO®	Ronde / blank	20 mm	Art.-Nr. / item-no.: D1022003
THERMEO®	Ronde / blank	25 mm	Art.-Nr. / item-no.: D1022006
THERMEO® SO	Ronde / blank	20 mm, A1	Art.-Nr. / item-no.: D1010125
THERMEO® SO	Ronde / blank	20 mm, A2	Art.-Nr. / item-no.: D1010126
THERMEO® SO	Ronde / blank	20 mm, A3	Art.-Nr. / item-no.: D1010127
THERMEO® MSI	Ronde / blank	16 mm, ozean-blau/ocean-blue	Art.-Nr. / item-no.: D1023002
THERMEO® MSI	Ronde / blank	20 mm, ozean-blau/ocean-blue	Art.-Nr. / item-no.: D1023003

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THERMEO® splint Cleaning and Care

Comprehensive documentation and important tips on the use and care of THERMEO® splints help you and the patients use THERMEO® products for a long and sustainable time.

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- dimensionally stable and accurate fit
- thixotropic, stable and controllable
- low shrinkage
- residue-free combustion

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LCR-1 is a light curing gel which can be used in a variety of applications. These applications include but are not limited to model cast technology, use as block-out gel for bleaching trays and as a repair material of 3D printed cast models. The gel is thixotropic, therefore steady and manageable when modeling. The marginal shrinkage of **LCR-1** leads to exceptional high dimensional stability and accuracy of fit.

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1 l Flasche / 1 l bottle

Art.-Nr. / item-no.: **U1900001**

New, water soluble cleaning agent for the removal of resin residues on additively manufactured dental molded parts.

- non explosive
- ideal for ultrasonic systems
- highly efficient (up to 5 times more than Isopropanole)
- ultrafiltrable
- compatible with **CLD-1**

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Exclusion of liability: We cannot be held responsible for false statements in this brochure.

Delivery of our products is effected according to our general terms of delivery and payment. Technical modifications reserved.

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pro3dure medical GmbH
pro3dure medical LLC