Natural beauty restored in one appointment



Lithium Disilicate CAD/CAM Block for chairside solutions





Since 1921 00 years of Quality in Dental

Natural beauty restored in one appointment

GC Initial[®] LiSi Block: new lithium disilicate block for one appointment dentistry


GC Initial[®] LiSi Block is a **fully crystallized lithium disilicate block** that delivers optimal physical properties without firing. This unique block features GC's proprietary **HDM** (High Density Micronization) **technology for CAD/CAM dentistry** to deliver high wear resistance, smooth margins and esthetic final results. This makes it an ideal, time saving solution for single visit chairside treatments.



- Save time, as no firing required
- Fully crystallized lithium disilicate
- Durable esthetic & accurate margins
- Natural opalescence


Just Mill, Polish and Place

GC Initial[®] LiSi Block can dramatically reduce process time: no need to fire, glaze, characterize and cool. This saves up to **40% in the time*** required to create your restorations, also reducing the chair time for you and your patient. You just need to mill, polish and place!



Source: GC R&D, Japan, Data on file. *Under testing conditions based on IFU.

Durable esthetics and smooth margins




Source: GC R&D, Japan, Data on file

Source: GC R&D, Japan, Data on file


- Source: GC R&D, Japan, Data on file
- Optimized acid and wear resistance to help preserve the esthetics of your restorations over time.
- Excellent edge stability for smooth margins.

More accurate margins

Being fully crystallized before milling, GC Initial[®] LiSi Block can be milled with smooth and accurate margins directly. Alternatively, it can be fired after staining and maintain great marginal accuracy.



GC Initial® LiSi Block GC Initial® LiSi Block (Polished) (Stain & glaze fired)



IPS e. max® CAD

Source: GC R&D, Japan, Data on file

Natural opalescence

GC Initial[®] LiSi Block is available in high translucency (HT) and low translucency (LT) and offers a natural opalescence in any light.

Choose your preferred finishing procedure


Superior gloss can be obtained in few minutes by polishing only, and the restoration is then ready for luting. For sophisticated aesthetic cases, remarkable results can be achieved with GC Initial® IQ Lustre Pastes NF and GC Initial[®] Spectrum Stains.



GC Initial[®] LiSi Block restoration under direct and


Courtesy of Dr. Javier Tapia Guadix, Spain





Inlays / Onlays

- Cavity wall angle: 6° with long axis
- Shoulder preparation



Full crowns

- Wall angle: 6~10°taper
- Deep chamfer or round chamfer preparation

Cement recommendation

Adhesive luting is recommended for GC Initial[®] LiSi Block. Both G-CEM ONE[™] and G-CEM LinkForce[®] from GC can be used for any type of indications using GC Initial[®] LiSi Block.



Function meets esthetics

«This case milled beautifully without marginal chipping and saved me a lot of chair time, since it doesn't require any firing! Just polished and cemented.»*

Karyn M. Halpern, DMD, MS Port Jefferson Station, NY



«GC Initial[®] LiSi Block has all of the advantages of lithium disilicate in terms of strength, esthetic, and bondability, without the need for firing.»^{*}

Yao-Lin Tang, DDS, San Mateo, CA





Courtesy of Karyn M. Halpern DMD, MS







Courtesy of Yao-Lin Tang, DDS






HDM technology for CAD/CAM dentistry



In 2016, with GC Initial[®] LiSi Press, GC introduced HDM (High Density Micronization) technology, which uses equally dispersed lithium disilicate micro-crystals to fill the entire glass matrix rather than using traditional larger size crystals. The clinical effectiveness of this technology has been proven after 5 years of clinical service¹⁾.

To bring fast solutions for one appointment dentistry, GC has further developed HDM technology for CAD/CAM dentistry by optimizing the crystal size and glass matrix stiffness. Thanks to this new technology, good machinability, marginal integrity, polishability, and wear resistance are achieved at the same time. The result is a strong and easy-to-mill block that offers the same strength with or without firing.



Source: GC R&D, Japan, Data on file.


Improved glass matrix stiffness for high mechanical

and high wear resistance

strength

Conventional lithium disilicate (IPS e. max[®] CAD[†])

HDM technology for CAD/CAM (GC Initial[®] LiSi Block)



Source: GC R&D, Japan, Data on file.


Workflow



Prepare



Scan



Design





Mill



Polish or characterize

Condition

Cement

Final result

Courtesy of Dr. Mark Kleive




Smaller crystal for easy milling



Ordering information

	GC Initial® LiSi Block CEREC®† Mandrel, size 14	
R man	Ref.	Shade
	012927	A1 HT - CEREC®†
Initial LiSi Block	012928	A2 HT - CEREC®†
Electron	012929	A3 HT - CEREC®†
	012930	B1 HT - CEREC®†
A A A A A	012931	A1 LT - CEREC®†
2.2.2	012932	A2 LT - CEREC®†
	012933	A3 LT - CEREC®†
	012934	B1 LT - CEREC ^{®†}



1) Cagidiaco EF, Sorrentino R, Pontoriero D, Ferrari M. 2020. A randomized controlled clinical trial on two types of lithium disilicate partial crowns. Am J Dent. 33(6):291-295.

Related products Image: Section ONE™ G-CEM ONE™ G-CEM LinkForce®

SKU #606035.09.2021.Ver3 • 800.323.7063 • www.gcamerica.com • www.gcamerica.com/training • © 2021 GC America Inc. ¹IPS e.max[®] CAD is a registered trademark of local Vivadent AG. CERC[®] and CERC[®] Tessera are registered trademarks of Dentsply Sirona. The dentists showcased in this brochure are compensated as Key Opinion Leaders for GC America. They have not been paid specifically to endows this product or provide these testimonials.

,'GC,'