Materials are our **DNA!**

pro **3d** ure
medical

D E N T A L

2 0 2 2

printo[®] dent

GR-22 flex

stay flexible

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"It has long been my wish to be able to produce flexible splints without impressions. With **printo[®] GR-22 flex**, we can now print our flexible splints in high quality. A significant step forward on my digital path."

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printo[®] GR-22 flex – the innovation for 3D printing.

The number of treatments of occlusal disorders by means of splint therapy is constantly increasing and thus the importance of the "dental splint" for a wide variety of indications such as bruxism, bite elevations or mouth protection. The demands placed on splint materials by patients, practitioners and dental technicians are manifold and range from durability, esthetics, wearing comfort and biocompatibility to processability and usability in digital processes. Based on biometric data, 3D printing processes can nowadays be used to produce complexly designed therapeutic splints efficiently. The latest member of the printo[®] family was developed with these aspects in mind – the innovative, flexible splint material printo[®] GR-22 flex with thermomemory effect, optimized for additive manufacturing. After all, who doesn't know it, the challenge of a "crooked impression" during conventional impression taking or imperfections during intraoral scanning that lead to "data holes" on the surface of a generated data set? A conventional, hard splint then wobbles on the plaster model and in the patient's mouth. The splint is subjected to greater stress at precisely these points. The risk of chipping or breaking is significantly increased. In addition, the wearing comfort for the patient is reduced. The innovative printo[®] GR-22 flex compensates for such inaccuracies thanks to its elastic properties combined with a thermomemory effect. printo[®] GR-22 flex becomes flexible at body temperature and thus adapts optimally to the tooth surface. This increases wearing comfort and minimizes the risk of fracture. printo[®] GR-22 flex thus supports the satisfaction of dental technicians, practitioners and patients through higher safety, exceptional wearing comfort and fewer "recalls".

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printodent® GR-22 flex is the flexible 3D printing resin for splints.

printodent® GR-22 flex gives you and your customers the flexibility you need! The development of the printodent® GR-22 flex was based on extensive clinical and material data from conventional, flexible splint materials. As a result, the printodent® GR-22 flex achieves increased flexibility in combination with a more pronounced thermomemory effect compared to its market competitors (Fig. 16), which is immediately noticeable when the splint is inserted in the mouth due to a higher pressure drop at the tooth. The wearing comfort for the patient is thus significantly increased. Would you also like to remain flexible in the selection and design of your "digital workflow"? A variety of validated processes based on a wide range of CAD softwares, 3D printers and post-treatment units are available for the printodent® GR-22 flex.

inside **exocad**

be outstanding

printodent® GR-22 flex – a resin with outstanding properties.

printodent® GR 22 flex meets the highest medical technology standards. Dental splints made of printodent® GR 22 flex have outstanding mechanical properties compared to their market competitors. For example, the standard specifications for impact strength (ASTM_D256) and elongation at break (ISO 37) are not only achieved, but exceeded by more than 3.5 and 2.5 times respectively (Fig. 17 and Fig. 18). The new development printodent® GR 22 flex is based on an innovative polymer matrix which ensures that the solubility of polymer components (Fig. 19) is also significantly lower than that of the comparative products. This ensures outstanding biocompatibility in this class of materials.

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printodent® GR-22 flex is durable and biocompatible.

Characterized by its exceptional wearing comfort, printodent® GR-22 flex has been developed for durability and reliability. The stability of the material was proven in an extensive chewing simulation study conducted by LMU Munich¹. Although printodent® GR-22 flex has a much higher flexibility than its market competitor, no significant differences in terms of material loss (volumetric and vertical) could be found between them (Fig. 20 and Fig. 21). printodent® GR-22 flex thus increases patient comfort without reducing clinical benefit! Furthermore, tested reworking and cleaning recommendations are available to sustainably ensure the special properties of a printodent® GR-22 flex splint over its entire service life.

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Fig. 20: printodent® GR-22 flex (orange) and KeySplint Soft (blue) volumetric degradation after 120,000 cycles chewing simulation, mm³.¹

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Fig. 21: printodent® GR-22 flex (orange) and KeySplint Soft (blue) vertical degradation after 120,000 cycles chewing simulation, mm.¹

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printodent® GR-22 flex is the innovative resin for 3D printing of splints with unique impact strength combined with high elongation at break. Its special viscoplastic properties in combination with a thermomemory effect combine highest comfort with highest biocompatibility.

- high elongation at break
- unique impact strength
- highest wearing comfort (usability studies pro3dure)
- lowest solubility/high biocompatibility
- thermomemory effect
- highest impact strength in its class
- "snap on" splints available

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pro3dure medical GmbH
Am Burgberg 13
58642 Iserlohn
Phone: +49 (0)2374 920050-0

USA

pro3dure medical LLC
Valley View Business Center
9825 Valley View Road
Eden Prairie, MN 55344
Phone: +1 (952) 426 1928

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Exclusion of liability: We cannot be held responsible for false statements in this brochure.

Delivery of our products is effected according to our general terms of delivery and payment. Technical modifications reserved.