PERMALINER

Denture Reliner

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
Fields of application:

Suitable for all PMMA resins.

- ▶ For fabrication of permanently soft linings in the dental laboratory.
- ▶ Removing pressure sores in the maxilla and mandible, adhesive difficulties and flabby ridges.
- ▶ For damming the A-line, as elastic margin for maxillary dentures and to secure the position of dentures through increased adhesion.
- ▶ For elastic covering of sharp-edged alveolar processes.
- ▶ For covering implant regions and as overdentures.
- ▶ For covering larger defects due to malformation or surgery.

Relining of new dentures with PERMALINER

Plaster cast

- ▶ First pour a plaster cast Class 3.

Spacer

- ▶ The spacer is best fabricated using a deep-drawing process. The thickness of the spacer is determined by the dentist and should be between 1 - 3 mm. The spacer can also be made of cold curing resin, i.e. tray material, and trimmed to the required thickness.

Investment

- ▶ Invest the wax set-up with spacer in a flask for press technology. Insulation of the teeth at the flask join is recommended with **FLEXISTONE®** or **FLEXISTONE® Plus**. Place flask in boiling water for approx. 4 min., then open and remove all traces of wax carefully. Remove spacer first.

Insulation

- ▶ Apply conventional alginate isolation. Then return spacer to its position.

Pressing and pre-polymerisation of the acrylate

- ▶ The acrylic should be mixed according to instructions and packed. A polyethylene foil is placed between the acrylic

and the spacer for the press process. Pre-press the flask at intervals of 10 - 15 min. at 1.5 tonnes. After opening the flask, remove the polyethylene foil and excess acrylic. Close flask with spacer and polyethylene foil and press again with clamp. Screw clamp tight, place flask in cold water and bring it to the boil. Boil for approx. 30 min., then allow to cool down. This pre-polymerisation prevents possible reactions with the acrylic monomer. Remove flask from the clamp and open. Remove spacer and polyethylene foil. Prior to packing **PERMALINER**, the acrylic must have attained a high degree of stiffness to avoid deformation.

Packing of PERMALINER

► Remove the ready-for-use **PERMALINER** from the tin with a clean spatula, apply to the acrylic and press slightly. Place a new polyethylene foil between the **PERMALINER** and the flask counter. Close the flask and pre-press at intervals at approx. 1.5 tonnes. After opening the flask, remove the polyethylene foil and excess **PERMALINER**. Reclose the flask and perform final pressing at approx. 3 tonnes. To polymerise the flask, place in clamp under pressure.

Polymerisation

► Place the clamped flask in a cold water bath and slowly heat to 100 °C.

Polymerisation time: boil for **1 hr. at 100 °C**.

Then allow flask to cool down slowly. Do not quench with cold water.

Polymerisation in the microwave unit

► Polymerisation is conducted in suitable commercially available microwave flasks (plastic, suitable for pressing) at 850 Watt within 10 min. Use dental plasters suitable for microwaving for investment. Investment and polymerisation must be conducted within 8 hours. Ovens without a rotating plate may require a gradual, all-side polymerisation by turning the flask. Deflask after self-cooling to room temperature.

Trimming of the denture

► Use normal procedure for trimming the acrylic. **PERMALINER** can be trimmed using cutters or grinding sleeves (15-20 thousand rpm) without heat build-up. Use pre-polishing discs for smoothing of rough spots.

Polishing of the acrylic

► Pre-polish and mirror polish acrylic employing standard methods. **PERMALINER** cannot be mirror polished. Lustrol gloss varnish is recommended as an alternative. The mixing ratio is 1:1. Prior to the application the base and functional margins have to be absolutely clean and dry. Apply **Lustrol gloss varnish** with a brush and allow to dry for approx. 10-15 min. **Lustrol gloss varnish** may only be applied to **PERMALINER**.

Relining of an existing acrylic denture with PERMALINER

► For relining a used denture with **PERMALINER**, the dentist takes a new impression of the existing denture base. Medium or low viscous impression materials are used for impression taking. The plaster cast, investment and insulation are carried out as already described.

The acrylic base

► Take the denture from the flask join and remove the impression material from the denture. Reduce the base by the desired Permaplast thickness (1.0 - 3.0 mm) by milling. To provide a smooth transition from acrylic to **PERMALINER**, it is recommended to cut an even circular groove into the acrylic in the transition zone to the desired Permaplast thickness. Remove any existing soft relining material completely.

► **Packing, pressing, polymerisation** and trimming of **PERMALINER** are carried out as already described.

Recommendations to avoid incorrect processing

1. Blistering:

- Insufficient pressure during pressing.
- Flask placed in the hot water bath.

- ▶ Packing of **PERMALINER** before the acrylic is stiff enough.
- ▶ Wetting of the acrylic with monomer liquid before applying **PERMALINER**.

2. Insufficient or no bonding:

- ▶ Base of denture not roughened deeply enough.
- ▶ The denture was wetted with alginate insulation.
- ▶ Hands and denture were not grease-free.
- ▶ Excessive heat build-up during trimming can lead to separation of **PERMALINER** at the margins.

Recommendations for repairs to **PERMALINER**

▶ **PERMALINER** can be repaired and supplemented by using Permaplast. To this purpose invest the denture. Roughen existing **PERMALINER** and resin well. Place fresh **PERMALINER** in the flask to supplement the repair zone and press (trial pressing with foil is possible). Then polymerise according to instructions for use.

Note

► **PERMALINER** may not be used in case of known allergies to any of the ingredients (see safety data sheet). Adverse events are not to be expected if used as instructed. Immune reactions, i.e. allergies, irritations, can however not be excluded as a matter of principle. If in doubt, it is recommended to perform an allergy test prior to use. Only to be used by trained specialist personnel.

DETAX is not liable for any damage caused by improper application.

► **Lustrol** gloss varnish: avoid contact of catalyst with eyes and skin, this can lead to irritation. In case of accidental contact, rinse eyes immediately with running water for at least 15 minutes, consult an ophthalmologist.

► Please follow the instructions of the safety data sheet!

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GmbH & Co. KG