



OPTIGLAZE color

EN	- Light-cured characterization material	2
DE	- Lichthärtendes Charakterisierungsmaterial	8
IT	- Materiale fotopolimerizzabile per caratterizzazione	14
FR	- Matériau de caractérisation photopolymérisable	20
ES	- Material para caracterización fotopolimerizable	26
PL	- Światłoutwardzalny materiał do charakteryzacji	32



EN - Characterization - Adjustment of Chroma



▲
Left : A3 + A-plus = A3.5
Middle : Original A3
Right : A3 + White = A2
(Ivory white)



▲
Left : B3 + B-plus = B4
Right : Original B3



▲
Left : C3 + C-plus = C4
Right : Original C3

EN - Characterization of Posterior Tooth


Occlusal side



Buccal side



Grey : Lavender = 5:1 (transparency)




Red brown

Pink orange

Red brown

Grey



Ivory white

Olive or Red brown (root)


EN - Characterization of Anterior Tooth



Ivory white : Pink orange
= 3 : 1
(crack liner)

Ivory white
(white band)

Ivory : Pink orange
= 4 : 1
(mamelon)



Simple technique



EN - Shades

A-Plus, B-Plus, C-Plus
For adjustment of shade (Chroma)



A-plus



B-plus



C-plus

For characterization of the fissures or reproduction of
mamelon, cervical area, white band, crack lines etc.



Yellow



Orange



Pink orange



White

(white band, chroma adjustment)



Ivory white



Pink



Red brown
(for fissures and proximals)



Olive
(for proximal areas)



Lavender
(for transparency)



Grey
(for transparency)



Blue
(for transparency)



Red
(characterizing of gum)

For diluting other shades and glazing (surface coating).



Clear




Clear HV (High Viscosity)

EN - Surface Characterization


Chairside

1. Preparation




- ① Finish with diamond bur or coarse silicone point. Oxygen inhibition layer must be removed.
- ② Rinse and dry.

2. Priming




When applying OPTIGLAZE color on hybrid ceramics (e.g. CERASMART) and composites (e.g. GRADIA/G-aenial), treat the surface using CERAMIC PRIMER II.

3. Characterize



Shake the bottle well and dispense in a dispensing dish, apply thinly on the prosthesis surface with a brush. Do not air blow.

4. Light Cure




See table 1. When using a gun-type light curing device, place the light tip as close as possible.

Table 1 : Light Curing Time

GC LABOLIGHT LV-III, II		5 min.
GC STEPLIGHT SL-I	Distance from light source <3cm	10 sec.
	Distance from light source >3cm	20 sec.
Halogen Light		40 sec.
Plasma Arc		8 sec.
LED Light (wavelength 400nm – 430nm)		40 sec.


Lab side

1. Preparation




- ① Sandblast with 25-50µm alumina (0.15MPa).
- ② Clean (steam/ultrasonic cleaner) and dry.

2. Priming




When applying OPTIGLAZE color on hybrid ceramics (e.g. CERASMART) and composites (e.g. GRADIA/G-aenial), treat the surface using CERAMIC PRIMER II.

3. Characterize



Shake the bottle well and dispense in a dispensing dish, apply thinly on the prosthesis surface with a brush. Do not air blow.

4. Light Cure



Set a 2cm platform in the curing device and place the prosthesis on it with the coated surface facing the light. Reverse the prosthesis and repeat the light curing to fully cure the undercut or areas in shadow.

EN - Internal Characterization

1. Preparation



Build up with Opaque and Dentin shade of high strength composite and light cure.

Note:


Do not remove the oxygen inhibition layer. Apply OPTIGLAZE color on the oxygen inhibition layer. If the oxygen inhibition layer is removed, treat the surface with CERAMIC PRIMER II and COMPOSITE PRIMER before applying OPTIGLAZE color.

2. Characterize



Shake the bottle well, dispense in a dispensing dish and apply to the prosthesis with a brush. Do not air blow.

3. Light Cure




Set a 2 cm platform in the curing device and place the prosthesis on it with the coated surface facing the light. Reverse the prosthesis and repeat the light curing to fully cure the undercut or areas in shadow.

Table 1 : Light Curing Time


GC LABOLIGHT LV-III, II	5 min.	
GC STEPLIGHT SL-I	Distance from light source <3cm	10 sec.
	Distance from light source >3cm	20 sec.
Halogen Light	40 sec.	
Plasma Arc	8 sec.	
LED Light (wavelength 400nm - 430nm)	40 sec.	

4. Build up with enamel



Continue building up with enamel shade on the oxygen inhibition layer.

5. Finishing & Polishing



Finish and polish using standard technique or apply OPTIGLAZE color CLEAR or CLEAR HV on the entire surface.

**GC EUROPE N.V.**

Researchpark
Haasrode-Leuven 1240
Interleuvenlaan 33
B-3001 Leuven
Tel. +32.16.74.10.00
Fax. +32.16.40.02.14
info@gceurope.com
www.gceurope.com

GC AMERICA INC.

3737 West 127th Street
USA-Alsip, ILL. 60803
Tel. +1.708.597.0900
Fax. +1.708.3771.5103
sales@gcamerica.com
www.gcamerica.com

GC GERMANY GmbH

Seifgrundstrasse 2
D-61348 Bad Homburg
Tel. +49.61.72.99.59.60
Fax. +49.61.72.99.59.66.6
info@germany.gceurope.com
www.germany.gceurope.com

GC ITALIA S.r.l.

Via Calabria
I-20098 San Giuliano Milanese
Tel. +39.02.98.28.20.68
Fax. +39.02.98.28.21.00
info@italy.gceurope.com
www.italy.gceurope.com

GC UNITED KINGDOM Ltd.

12-15, Coopers Court
Newport Pagnell
UK-Bucks. MK16 8JS
Tel. +44.1908.218.999
Fax. +44.1908.218.900
info@uk.gceurope.com
www.uk.gceurope.com

GC FRANCE s.a.s.

9 bis, Avenue du Bouton d'Or – BP 166
F-94384 Bonneuil sur Marne Cedex
Tel. +33.1.49.80.37.91
Fax. +33.1.49.80.37.90
info@france.gceurope.com
www.france.gceurope.com

GC IBÉRICA

Dental Products, S.L.
Edificio Codesa 2
Playa de las Américas, 2, 1º, Of. 4
ES-28290 Las Rozas, Madrid
Tel. +34.916.364.340
Fax. +34.916.364.341
info@spain.gceurope.com
www.spain.gceurope.com

GC AUSTRIA GmbH

Tallak 124
A-8103 Rein bei Graz
Tel. +43.312.45.40.20
Fax. +43.312.45.40.20.40
info@austria.gceurope.com
www.austria.gceurope.com

GC BENELUX B.V.

Edisonbaan 12
NL-3439 MN Nieuwegein
Tel. +31.30.630.85.00
Fax. +31.30.605.59.86
info@benelux.gceurope.com
www.benelux.gceurope.com

GC EUROPE N.V.

East European Office
Siget 19b
HR-10020 Zagreb
Tel. +385.1.46.78.474
Fax. +385.1.46.78.473
info@eeo.gceurope.com
www.eeo.gceurope.com

GC NORDIC AB

Finnish Branch
Vanha Hommaksentie 11B
FIN-02430 Masala
Tel. & Fax. +358.9.221.82.59
info@finland.gceurope.com
www.finland.gceurope.com

GC NORDIC AB

Danish Branch
Harbour House
Sundkrogsgade 21
DK-2100 København
Tel. +45 23 26 03 82
info@denmark.gceurope.com
http://nordic.gceurope.com

GC NORDIC AB

Box 70396
107 24 Stockholm
Sweden
Tel. +46.8.506.361.850
info@nordic.gceurope.com
http://nordic.gceurope.com

GC AUSTRIA GmbH

Swiss Office
Bergstrasse 31
CH-8890 Flums
Tel. +41.81.734.02.70
Fax. +41.81.734.02.71
info@switzerland.gceurope.com
www.switzerland.gceurope.com