LC-GEL

INSTRUCTIONS

DESCRIPTION GÉNÉRALE

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GENERAL DESCRIPTION

LC-Gel is a light-curing material ready to use for modelling of inlays, onlays, veneers, secondary parts, other modelling constructions and for blocking out undercuts. LC-Gel is applied directly from the syringe.

COMPOSITION

LC-Gel is a modelling gel MMA-free based on multifunctional acrylates. Gel cured with a standard light-curing unit of 320-500 nm. LC-Gel burns out free of residues at 480 °C (900 °F).

METHOD OF USE

Preparation of the model

The model must be free of dust and grease. It is necessary to isolate the model with a separating liquid, if the part has to be taken out from the model.

Application

Remove the cap of the LC-Gel syringe and put a dispensing tip. Apply LC-Gel on the desired area to a maximum thickness of 3 mm and cure with a curing-light of 320-500 nm for 5 minutes.

Note : LC-Gel polymerizes under the influence of daylight. Do not extrude the gel from the syringe unless you intend to use it immediately.

Too short exposure time or low intensity of light may lead to a loss of stability.

STORAGE

To be stored below 25 °C. Avoid direct sunlight. Excess material can be used again if stored under light protection.

To Order: DENPLUS Inc.
333-M Chemin du Tremblay, Boucherville, QC, Canada, J4B 7M1
Tel.: (450) 641-1330 or toll-free 1-888-344-4424
Fax: (450) 641-1337
www.denplus.com

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Pour commander : DENPLUS Inc.