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

LED light curing light

Introduction

This LED light can be used to polymerize light-curing resins between 420 nm and 480 nm. The lithium battery can be used between 80 to 100 times at full power between each recharge.

Installation

IMPORTANT: To prevent the lamp from turning on during transport. A sticker has been affixed to the battery. Unscrew the base of the lamp and remove the small colored sticker on top of the battery. Insert the lamp head containing the optical fiber. Attention, optical fiber is fragile. Install the orange shield.

Plug the charger into the charging unit and place the lamp in the charging unit.

Specification

Dimension: 24 x 200 mm

Weight: 220 g

Voltage: 110V - 220VAC, 50-60 Hz

Battery: Lithium, 2200 mAh

Operation

Switch on the lamp by pressing the button P. Close the lamp by pressing the P button again. The lamp will switch off automatically after one minute if no action is taken.

The M button allows you to choose between 3 modes: continuous, gradual and intermittent. The continuous mode gives full power for the desired time. The gradual mode increases the power gradually for about 5 seconds and maintains full power for the end of the cycle. The intermittent mode closes and opens the light regularly during the desired time. Every 5 seconds, the lamp emits a beep.

To change the time, hold the M button in the selected mode until the indicator light flashes. Change the time in 5-second intervals by pressing the P button. The maximum time is 40 seconds for each operating mode.

Sterilization of the optical fiber: Soak a few minutes in a 75% alcohol solution.