## Intended use

EN

Material on methacrylate resin basis for DLP systems with 385nm resp. 405nm-LED for manufacturing of dental surgical drill guides for tooth-supported and/or teethridge-supported and surgical invasive application.

**Indications for use:** Manufacturing dental drill guides for placing dental implants.

**Indication:** Generative manufacturing of drill guides using 3D printing. This serves the dentist as a guide rail during drilling.

**Contraindication:** Only apply FotoDent® guide intraorally in a fully polymerized state. If you are allergic to acrylic, FotoDent® guide should not be used.

**Patient target group:** Patients with a single-tooth gap, partially and/or completely edentulous upper and/or lower jaw.

**Intended users:** The dental technician or specialist staff in the dental laboratory is responsible for the production of the drill guide in 3D printing. The dentist uses the finished drill guide on the patient.

## Processing

- Alignment of drill guide: vertical, supports, connection of supports to a baseplate.
- During the application of FotoDent® guide we recommend wearing personal protective equipment e. g. suitable gloves, safety goggles etc.
- Please see further steps within printer manual.
- Carefully pour FotoDent® guide into the pre-determined container of the production unit.
- Remove all bubbles with a cleaned object.

Processing temperature	
Rapid Shape D and DII series	23 ± 3 °C
Asiga MAX UV	35 ± 3 °C

- You can find more qualified printers on our website [www.dentamid.dreve.de](http://www.dentamid.dreve.de).
- Select the buildstyle for FotoDent® guide in the printer software. Make sure that all software is up to date.
- After completion of the building process, direct further processing (see following steps) is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.

- Cleaning and post curing incl. supports and baseplate to avoid warpage.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 2x6 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended.
- For post-curing of the parts, please select the parameter settings for your process.

Post-curing 385 nm	
PCU LED N <sub>2</sub>	
Rapid Shape DII / Asiga MAX UV	5 minutes, 33 % light output (incl. nitrogen)
FotoDent® flash	
Rapid Shape DII / Asiga MAX UV	6 minutes + 1 minute gas pre-flow
RS Cure (Gas)	According to respective parameter set for FotoDent® guide 385 nm.
Post-curing 405 nm	
FotoDent® flash	
Rapid Shape D	10 minutes + 1 minute gas pre-flow

- After post-curing, the supporting structures and base plate must be removed.
- The surface can be smoothed by sanding or polishing.
- Fitting of suitable drill sleeves (from other manufacturers) into the drill guide. Please refer to the instruction for use of the supplier/manufacturer for the adjustment of the drill sleeves.
- Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate FotoDent® guide. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1–2 hours with FotoDent® guide.

**To avoid damage of material quality do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may impair biocompatibility, lead to modified mechanical characteristics and/or color variations of the FotoDent® guide material.**

## Possible dangers

Please note the safety instructions for processing FotoDent® guide.

**Hazard statements: H319** Causes serious eye irritation. **H317** May cause an allergic skin reaction. **H412** Harmful to aquatic life with long lasting effects. (REF D35600 only) **H413** May cause long lasting harmful effects to aquatic life. (REF D35650 only)

**Precautionary statements: P261** Avoid breathing dust/fume/gas/mist/vapours/spray. **P273** Avoid release to the environment. **P280** Wear protective gloves/eye protection/face protection. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P333+P313** If skin irritation or rash occurs: Get medical advice/attention. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

**Contains:** 2-hydroxyethyl methacrylate, hydroxypropylmethacrylate, monoester with propane-1,2-diol, Diphenyl(2,4,6-trimethylbenzoyl) phosphin oxide

## Cleaning and disinfection

For cleaning and disinfection, we recommend the use of ethanol (purity ≥ 60 %) or a 0.2 % chlorhexidine solution, each with a contact time of 5 minutes, then allow to dry well.

## Sterilisation

Do not use thermal methods for sterilisation, otherwise deformation of the surgical guide cannot be ruled out. FotoDent® guide cannot be sterilised.

## Attention

Select FotoDent® guide material that matches the wavelength of the used printer (385nm or 405nm).

## Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® guide packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

## Serious incidents

All serious incidents occurring in relation to the product shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is settled.

## Note for single use

Due to its physical features, the printed and cured FotoDent® guide can only be used for one drill guide.

## Disposal

Disposal of contents/container in accordance with local/regional/national/international regulations and according to the safety data sheet.

## Glossary

	Manufacturer		Protect from sunlight
	Batch Code		Temperature limitation
	Item number		Expiry date
	Notified body registration number		Use by only professionals
	Medical device		Please note instructions for use

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