


DB Plus COMPOSITE POSTERIOR TEETH
DB Plus COMPOSITE POSTÉRIEURES

20° COMPOSITE POSTERIOR TEETH
20° COMPOSITE POSTÉRIEURES



30° COMPOSITE POSTERIOR TEETH
30° COMPOSITE POSTÉRIEURES



db plus COMBINATION CHART

MAXILLARY ANTERIOR MOULD	ARTICULATES WITH MANDIBULAR MOULD	COMBINES WITH 0° POSTERIOR MOULD	COMBINES WITH F-10 POSTERIOR MOULD	COMBINES WITH 20° POSTERIOR MOULD	COMBINES WITH 30° POSTERIOR MOULD
O40	L1	29	F29	29M or 29L	S30/L30/31
O45	L1	29	F29	29M or 29L	S30/L30/31
O52	L1	29	F29	29M or 29L	S30/L30/31
O55	L8	29	F29	29M or 29L	S30/L30/31
O65	L9	31	F31	31M or 31L	S32/L32/33
O75	L13	31	F31	31M or 31L	S32/L32/33
O80	L13	33	F33	33M	33/X34
O89	L13	33	F33	33M	33/X34
O99	L16	33	F33	33M	33/X34
T40	L1	29	F29	29M or 29L	S30/L30/31
T43	L3	29	F29	29M or 29L	S30/L30/31
T53	L3	29	F29	29M or 29L	S30/L30/31
T54	L5	29	F29	29M or 29L	S30/L30/31
T63	L6	31	F31	31M or 31L	S32/L32/33
T76	L11	31	F31	31M or 31L	S32/L32/33
T81	L9	33	F33	33M	33/X34
T90	L15	33	F33	33M	33/X34
R40	L1	29	F29	29M or 29L	S30/L30/31
R41	L1	29	F29	29M or 29L	S30/L30/31
R58	L3	31	F31	31M or 31L	S32/L32/33
R61	L5	31	F31	31M or 31L	S32/L32/33
R69	L9	33	F33	33M	33/X34
R78	L15	33	F33	33M	33/X34
R88	L13	33	F33	33M	33/X34
R91	L15	33	F33	33M	33/X34
S45	L1	29	F29	29M or 29L	S30/L30/31
S56	L5	29	F29	29M or 29L	S30/L30/31
S65	L7	31	F31	31M or 31L	S32/L32/33
S67	L7	31	F31	31M or 31L	S32/L32/33
S71	L11	31	F31	31M or 31L	S32/L32/33
S79	L13	31	F31	31M or 31L	S32/L32/33
S83	L13	33	F33	33M	33/X34
S98	L13	33	F33	33M	33/X34

myerson.

MANUFACTURER
Myerson Company Limited
3 trinity avenue
Laventille, Trinidad & Tobago
phone 868 623 1007

DISTRIBUTED BY
Myerson LLC
5106 North Ravenswood
Chicago IL 60640-2713
USA
phone 800 423 2683

EC REP
EUROPEAN REPRESENTATIVE
Obelis s.a
Bd. Général Wahis 53
1030 Brussel, Belgium
phone +(32) 2.732.59.54
fax +(32) 2.732.60.03
email: mail@obelis.net

Preformed plastic denture teeth for use in making temporary or removable dentures. Documented as per ISO 22112:2017
2020-01-29

Make great removables
with Myerson teeth
& Myerson's advanced
thermoplastic materials

DURAFLEX
PARTIALS



VISICLEAR
CLASPS & FRAMEWORKS



DURACETAL
COSMETIC APPLIANCES



ema
from
myerson.

With 4 strap strengths and
9 strap lengths, a custom
EMA treats snoring and
obstructive sleep apnea
safely and effectively, giving
you the sleep you need to
lead a full, healthy life.



ema is a trademark of frantz design inc.
US Patent Numbers: 6,106,265; D611,153 & Patent Pending

myerson DB Plus
mould chart

myerson.

**We make the teeth
that make the smile.**

DB Plus RESIN ANTERIOR TEETH • MAXILLARY
DB Plus RÉSINE ANTÉRIEURES • DU MAXILLAIRE SUPÉRIEUR

OVOID EN FORME DE PELLE



TAPERING TRIANGULAIRE



DB Plus RESIN ANTERIOR TEETH • MAXILLARY
DB Plus RÉSINE ANTÉRIEURES • DU MAXILLAIRE SUPÉRIEUR

RECTANGULAR RECTANGULAIRE



SQUARE CARRÉ



DB Plus RESIN ANTERIOR TEETH • MANDIBULAR
DB Plus RÉSINE ANTÉRIEURES • DU MAXILLAIRE INFÉRIEUR



DB Plus COMPOSITE POSTERIOR TEETH
DB Plus COMPOSITE POSTÉRIEURES

0° COMPOSITE POSTERIOR TEETH • 0° COMPOSITE POSTÉRIEURES



F10 COMPOSITE POSTERIOR TEETH • F10 COMPOSITE POSTÉRIEURES



myerson.