



PRINTER VALIDATION PROCESS

CERTIFIED
VALIDATED
PROCESS
RELIABLE



CURRENT
VALIDATION
TABLE

WWW.DETAX.COM



FREEPRINT®
MED

FREEPRINT®
TEC

CO-BRANDING
MED & TEC

Qualification	385 nm														405 nm										
	ASIGA Max / Mini	ASIGA Pico2	ASIGA PRO2	ASIGA PRO 4K	Rapid Shape D10 / D20 Series	Rapid Shape D30 / D40 Series	Rapid Shape D50+ Series	Rapid Shape D70 / D90 Series	Straumann P Series	Ivoclar PrograPrint PRS	MICROLAY Versus	Miicraft Alpha Series	Miicraft Prime / Hyper Series	Miicraft Ultra Series	Miicraft Profession / Advance Series	W2P SoftFlex	Ackuretta SOL	MicroRay Eye Pro	Phrozen Sonic 4K / XL	Prusa Research MEDICAL ONE	Shining 3D Accu-Fab D1s	Shining 3D Accu-Fab L4D	Dentsply Sirona Primeprint	Shera SHERAprint Series	
TEMP	✓	✓	✓	✓	✓	✓	□	✓	✓	✓	✓	✓	✓	✓	✓	✓	□	●	□	●	●	●	●		
CROWN	✓	●	●	✓	✓	✓	□	✓	✓	✓	✓	□	□	□	□	□	□	●	□	●	●	●	●		
DENTURE	✓	✓	✓	✓	✓	✓	□	●	✓	✓	✓	✓	✓	✓	✓	✓	□	●	□	●	●	●	●		
DENTURE IMPACT	□	●	●	□	□	□	□	●	□	□	□	□	●	●	●	●	●	●	●	●	●	●	●		
ORTH0	✓	✓	✓	✓	✓	✓	□	✓	✓	✓	✓	✓	✓	✓	✓	✓	●	●	□	●	●	●	●		
SPLINT 2.0	✓	✓	✓	✓	✓	✓	□	●	✓	✓	✓	✓	✓	✓	✓	✓									
SPLINTMASTER	✓	●	●	✓	✓	✓	□	●	✓	✓	□	□	□	□	□										
IBT	✓	✓	✓	✓	✓	✓	□	●	✓	✓	✓	✓	✓	✓	✓	✓									
TRAY 2.0	✓	✓	✓	✓	✓	✓	□	●	✓	✓	✓	✓	✓	✓	✓	✓		●	●	●	●	●	●		
TRYIN	✓	✓	✓	✓	✓	✓	□	✓	✓	✓	✓	✓	✓	✓	✓	✓	●	●	□	●	●	●			
MODEL	✓	✓	✓	✓	✓	✓	□		✓		✓		✓	●	●	✓		✓	✓						
MODEL 2.0	✓	✓	✓	✓	✓	✓	□	●	✓		✓	□	✓	●	●	✓		✓	✓	□	✓	✓			
MODEL PRO	✓	●	●	✓	□	□	□	●	□		□	□	□	●	●			□	□	□	✓	✓			
MODEL T	✓	✓	✓	✓	✓	✓	□	✓	✓		✓	□	✓	●	●	✓		✓	✓	□	✓	✓			
MODEL WW	✓	✓	✓	✓																					
GINGIVA	✓	✓	✓	✓	✓	✓	□	●	✓	✓	✓	□	✓	●	●	✓		✓	✓	□	✓	✓			
CAST 2.0	✓	✓	✓	✓	✓	✓	□	●	✓	✓	✓	□	✓	●	●	✓		✓	✓	□	✓	✓			
Primeprint (made by DETAX)																							✓		
SHERAprint Splint Taff/Flex (made by DETAX)																								✓	