


FREEPRINT®

»» SPLINTMASTER «« TAFF & FLEX

Ultimate solution for splints, repositioners, mouth- & nightguards

NEW



- ✓ Extreme impact strength
- ✓ Flexible and robust
- ✓ Oral resistant
- ✓ Long-term colour stable
- ✓ Odour- & tasteless



FREEPRINT® SPLINTMASTER

Light-curing 3D resin for printing of flexible splints. Available in two levels of flexibility: TAFF & FLEX.

The clear-transparent material is characterized by extreme impact strength. Splints made of FREEPRINT® SPLINTMASTER are biocompatible, long-term stable and robust, without becoming brittle. Thanks to highest elastic recovery the splints always return to their original shape. Perfectly tailored for tension-free wear and easy insertion & removal. The clear-transparent splints are odour- and tasteless and free from MMA, THF-MA and TPO.

FREEPRINT® SPLINTMASTER meet all requirements for a **Class IIa Medical Device according to the new MDR**, has **U.S. FDA 510(k)** and **Health Canada Class II license**.

- ✓ Tension-free wear
- ✓ Effortless insertion & removal
- ✓ Easy to clean
- ✓ Clear-transparent
- ✓ MMA-, THF-MA- & TPO-free



REF 04433
REF 04432

FREEPRINT® SPLINTMASTER TAFF, 1.000g
FREEPRINT® SPLINTMASTER FLEX, 1.000g

