

A5601 Stump separating unit CUT



Product group

Machines for cutting operations of plaster or investment patterns

Main features

- ⇒ Diamond-coated separating disc
- ⇒ Pivoting arm with dual safety switch
- ⇒ Transparent protection shield on the blade
- ⇒ Adjustable lighting unit
- ⇒ Model-holder locking magnetic base
- ⇒ Model-holder with interchangeable plates according to the different working techniques
- ⇒ Electric feeding socket and single port for connection to the dust extraction system

Gearmotor

Gearing system with toothed driving belt

Speed of rotation

Fixed at 2800 rev/min

Direction of rotation

Clockwise

Blade dimensions

Ø 80mm – 0.25mm in thickness

Voltage line

230v AC 50/60 Hz

Absorption

100W

Width

230 mm

Depth

500 mm

Height

360 mm

Net weight

9.3 Kg

Gross weight

14.4 Kg