

A5801 VACUUM PUMP



Product group	Vacuum equipment
Main features	<ul style="list-style-type: none"> ⇒ Essential accessory to create vacuum in the heating chamber of porcelain furnace; it can serve also other equipment ⇒ High flow and power ⇒ Extremely silent ⇒ No maintenance required
Vacuum pump	Double piston "dry type" (oil free)
Vacuum level	745 mmHg
Flow	50 l/min
Connecting pipe	Ø 6x11 mm – lenght 1,6 m
Line voltage	230v A.C. 50 Hz
Absorption	200W
Width	150 mm
Depth	300 mm (including the space required for bending the pipe)
Height	200 mm
Net weight	7.6 kg
Gross weight	8.3 kg