

# New Design, Same Great Value



## FireLite Porcelain Furnace

- **NEW** Integrated Control Board Design
- **50 Firing Programs** customizable for any type of Porcelain
- **2-Point Factory Calibration** ensures firing accuracy for both low fusing and high fusing porcelains
- **Fast Cool Feature** uses the air from the muffle fed through a vacuum to reduce cycle time
- **Cooling Tray** accessory is large and elevated which allows for easy staging of work and quickly cools just fired restorations
- **Single Stage Furnace** crystallizes E-max with standard crystallization program. Speed crystallization requires a Two Stage furnace like the Pro200
- **Customizable Sound** to differentiate between multiple furnaces

**Item #91057** 100 VAC ± 10%, 50/60 Hz  
**Item #91050** 115 VAC ± 10%, 50/60 Hz  
**Item #91055** 230 VAC ± 10%, 50/60 Hz



**Whip Mix Corporation** 361 Farmington Avenue 800.626.5651  
 Louisville, KY 40209 502. 637.1451  
 www.whipmix.com Fax 502.634.7169



Connect With Us!