



GC FUJIROCK



FUJIROCK	FUJIROCK EP	FUJIROCK EP PREMIUM	NEW FUJIROCK	NEW FUJIROCK IMP	NEW FUJIROCK FAST SET
GYPSUM DIE, TYPE 4	All kinds of prosthetics	Implant-ceramic complex	All kinds of prosthetics	For implant CAD-CAM	For high precision prosthetics
SETTING EXPANSION	≤ 0.09%	≤ 0.09%	≤ 0.08%	≤ 0.06%	≤ 0.08%
COMPRESSIVE STRENGTH	≈ 53 MPa	≈ 60 MPa	≈ 58 MPa	≈ 54 MPa	≈ 58 MPa
SETTING TIME	10-13 min.	9-12 min.	11 min.	11 min.	5 min.
REMOVAL TIME	30-40 min.	25-35 min.			
	12 kg	11 kg 4 kg	3 kg	3 kg	3 kg
AVAILABLE COLORS	GB	PY, PW, IB, TG	I, GB	I, GB	GB

*All types of Fujirock are treated to be scannable

*Fujirock premium (**PW, IB, TG**) colors are available just in 4 kg size

***I** = IVORY, **G** = GOLDEN, **PY** = PASTEL YELLOW, **PW** = POLAR WHITE, **GB** = GOLDEN BROWN, **IB** = INCA BROWN, **TG** = TITANIUM GREY

*Prices can be changed without notice

www.denplus.com

1-888-344-4424