

Technical data sheet

These Phare Médica nitrile gloves are ideal for use in the areas of dental, cosmetic and medical examinations. Comfortable and resistant, it offers superior dexterity and optimal protection for users.

Generals informations

Trade name	Phare Médica nitrile glove for medical examinations
Brand	Phare Médica
Product code	777-320XS (Extra Small), 777-320S (Small), 777-320M (Medium), 777-320L (Large), 777-320XL (Extra Large)
Application	Disposable nitrile gloves for medical examination
Manufacturing country	China
Distributed by:	Phare Médica 1660 Boulevard Saint-Paul, Local 100, Saguenay, G7J 4N1, Québec, Canada

Specifications

Material	Stretch nitrile
Color	Turquoise
Dimension	Thickness : Palm 3.15 Mils Fingertips : 3.54 Mils Cuff : 2.76 Mils Thickness : Average 3.15 Mils Length : 9.5 inches
Weight	XS : 2.8 grams , S : 2.8 grams , M : 3.0 grams , L : 3.5 grams , XL : 4.0 grams
Model	Ambidextrous
Hypoallergenic	Latex free / Protein free
Other informations	Textured fingertips for a secure grip Powder free Non-sterile

Chemotherapy test result - D6978

Chemotherapy drug tested	Breakthrough detection time	Chemotherapy drug tested	Breakthrough detection time
Cisplatin (1mg/ml)	>240 Minutes	Fluorouracil (50mg/ml)	>240 Minutes
Dacarbazine (10mg/ml)	>240 Minutes	Pacilitaxel (6mg/ml)	>240 Minutes
Etoposide (20mg/ml)	>240 Minutes	Mitomycin-C (0.5mg/ml)	>240 Minutes
Cyclophosphamide(20mg/ml)	>240 Minutes	Vincristine Sulfate (1mg/ml)	>240 Minutes
Doxorubicin HCL (2mg/ml)	>240 Minutes	Thiotepa (10mg/ml)	Non recommandé
Carmustine (BCNU) (3.3mg/ml)	Non recommandé		

Qualification

ASTM D-6319

ASTM D-6978

ASTM F-1671

EN 455 et EN 374

Food Save

Technical data sheet

Packaging

Packaging	100 gloves / box, 10 boxes / case
Box measurements	23 cm X 13 cm X 6.5 cm
Box weight	600 grams

Others informations

Storage conditions	Store in a dry, cool and well ventilated place. Do not store above 37°C
Caution	The components used in the manufacture of gloves may cause an allergic reaction in some users. Follow your establishment's usage policies. People with known sensitivity should avoid contact.

