

A1211 FORCED VENTILATION FUME SUCTION SYSTEM (suitable for CF-2 BURN-OUT FURNACES)



Product group	Fume suction system for burn-out furnaces (it can not be used as self-operating unit)
Main features	<ul style="list-style-type: none"> ⇒ Indispensable accessory to pipe the fumes originated during the baking stage of the casting rings in case the Laboratory is not fitted with extraction chimney ⇒ It can be added-on to every size of CF-2 burn-out furnaces
Automatic operation	YES
Line voltage	230V A.C. 50 Hz (different voltage upon request)
Width	150 mm
Depth	170 mm
Height	600 mm
Exhaust pipe diameter	60 mm
Net weight	3,3 Kg
Gross weight	4,0 Kg