

A3502 "ISO"



Product group

ISO SURVEYOR/PARALLELOMETER complete

Main features

- ⇒ **Double-joint arm supporting any micromotor, turbine or hot spatula**
- ⇒ **Most suitable to carry out milling operations both on wax and or metal**
- ⇒ Complete with model-holder base
- ⇒ Stainless steel working plan
- ⇒ Complete set of accessories allowing to carry out:
 - Restrained vertical movement for drilling and boring
 - Restrained horizontal movement for coulisse
 - locking the articulated joints
- ⇒ Stainless steel tool storage compartment
- ⇒ At same time, all procedures usually carried out with the manual Surveyor can be easily executed with this instrument

Number of arms

1 - double-joint arm fitted with bivalent head, featuring on one side the micromotor locking pliers and on the other side the slide-guide self-lubricating bush for the pliers holding rod

Width

190 mm

Depth

270 mm

Height

330 mm

Net weight

5,2 Kg