

A3503C ISO Colou Top



Product group
Main features

ISO SURVEYOR/PARALLELOMETER complete

- ⇒ **Double-joint arm to support any micromotor, turbine or hot spatula**
- ⇒ **Most suitable to carry out milling operations both on wax and on metal**
- ⇒ Complete with model-holder base
- ⇒ Stainless steel working plan
- ⇒ Complete set of accessories allowing to carry out:
 - Restrained vertical movement for drilling and boring
 - Restrained horizontal movement for coulisse
 - locking the articulated joints
- ⇒ Stainless steel tool storage compartment
- ⇒ At same time, all procedures usually carried out with the manual Surveyor can be easily executed with this instrument
- ⇒ Built-in lighting and cooling systems
- ⇒ Special version featuring enameled base and satined or oxidized light alloy components available in blue, red or natural

Number of arms

1 - double-joint arm fitted with bivalent head, featuring on one side the micromotor locking pliers and on the other side the slide-guide self-lubricating bush for the pliers holding rod

Lighting

12V Revolving halogen lamp

Line voltage

230V A.C. 50Hz (different voltage upon request)

Absorption

6W

Width

200 mm

Depth

340 mm + distance to the electric cable and the air pipe

Height

330 mm

Net weight

7,0 Kg