

A5090 HIGH SPEED TURBINE ZIRKO-SURFIN 300k



Product group Laboratory air turbine with water spray

Main features ⇒ Powerful and silent, best suited for all finishing operations on æsthetic

products in ceramics, litium silicate or zirconia

 \Rightarrow The rotor at 2,5-3 BAR reaches 300.000 rpm

⇒ Refrigeration of both the tool and the machined surface is ensured by a constant flow of air and can be increased with the adition of water spray

⇒ Maintenance and lubrication are very simple

⇒ Control unit incorporating everything needed for air treatment, any adjustments and the tank for distilled water

⇒ Extremely small and lightweight handpiece

⇒ Circular footswitch, operated in any position

Rotating speed 300.000 rpm

Rotation Fixed
Foot control YES

Standard pliers Ø 1.6 mm

Tool needed to change the bur YES

Handpiece dimensions Ø 17 mm – lemgth 130 mm

Handpiece weight 85 g

Supply Air supply 2,5-3 BAR

Air consumption 50-60 l/min

Line tension N.A.

Refrigerating liquid Distilled water

Tank capacity 170 ml

Overall dimensions of supply unit

(WxDxH)

260x110x190 mm

Dimensions of foot control Ø 150 mm Lenght of pipes (foot control) 1.8 m

Weight 2,7 – 3,8 Kg (Gross weight)