

## A6003D HOT 2 Plus Waxing Unit



Product group

**Electrical appliances for hot waxing-up**

Main features

- ⇒ Extremely compact in size
- ⇒ Electronic temperature control (microprocessor)
- ⇒ Wide temperature range from 50°C to 200°C, allowing to choose the best temperature according to every type of wax product, preserving its properties.
- ⇒ 2 independent working places
- ⇒ 2 recordable temperatures for each working place
- ⇒ 3-digit-display showing both the set temperature and ultraflat keyboard
- ⇒ Conversion key (Celsius or Fahrenheit scale)
- ⇒ Thermally insulated handpiece with tool anti-rotation system
- ⇒ 2 handpiece bench holders
- ⇒ Tool seat directly on the control unit + practical rear compartment
- ⇒ Two available sets of tools: 6 modelling tips, and 4 milling tips

Line tension

230v A.C.

Handpiece tension

24v

Absorption

20W

Length

130 mm

Width

110 mm

Height

60 mm

Net weight

1,0 Kg (1,6 kg packed)