



## LC-GEL

Light-curing resin gel  
DenPlus.

Syringe of 3 g and 2 tips  
Color: **Red**



LC-GEL is a ready to use modelling resin, fabricated of multifunctional light curing acrylates. Its properties make it multi-usages and easy to use as well.

- Ready to use with a very stable consistency
- Wide range of use: Crowns & Bridges, Attachments, Implants, Metal soldering, Inlays & Onlays, Indirect Posts
- Cured with a standard light cure in few seconds
- Burns off without leaving any residues



**Substitute to GC Pattern Resin LS to avoid mixing problems since its ready to use gel with good consistency**

Distributed exclusively by DenPlus

[www.denplus.com](http://www.denplus.com)

1-888-344-4424