



# DP Estética PMMA discs

Powered by **UNIDESA**



DP Estética discs are the choice of dental professionals  
**who value product quality over all else!**

DP Estética PMMA discs are manufactured in partnership with Union Dental SA using only the highest quality raw materials and exacting quality control processes, resulting in outstanding hardness and wear resistance, shade precision and aesthetics.

## ***All Estética PMMA discs:***

- Use highly cross-linked and reinforced PMMA with high molecular weight resulting in superior hardness and wear-resistance
- Use 100% mineral pigments resulting in excellent colour stability and repeatable shade accuracy
- Meet rigorous quality control ISO standards. 100% of DP Estética PMMA discs are inspected off the line and each disc is manufactured individually
- Are easy to mill and polish due to the quality and density of the raw material



unidesa · *odi*



**1-888-344-4424**  
[www.denplus.com](http://www.denplus.com)



DP Estética discs are the choice of dental professionals  
**who value product quality over all else!**



### DP Estética 10-Layer tooth coloured PMMA discs

When ultimate aesthetics are the priority.

- Milled dentures
- Long-term aesthetic provisionals
- Perfect gradation of colour
- Similar fluorescence to natural teeth

**98.5mm | Shades:** A1, A2, A3, B1, C2, BL1  
20mm: \$112.75



### DP Estética 8-layer translucent tooth coloured PMMA discs

Unique formulation with 40% greater translucency than Estética 10L.

- Milled dentures, long-term aesthetic provisionals
- Perfect gradation of colour
- Similar fluorescence to natural teeth

**98.5mm | 16mm Shades:** A1, A2, A3, A3.5, A4, B1, C1, C2, C3, C4, D2, D3, D4, BL1 - \$84.50  
**98.5mm | 20mm Shades:** A1, A2, A3, B1, C2, BL1  
\$96.25

### DP Estética 5-layer tooth coloured PMMA discs

Ideal for short-term temporaries or mock-ups where strength, aesthetics and shade precision are valued.

**98.5mm | Shades:** A1, A2, A3, B1, C1, D2  
16mm: \$84.50 | 20mm: \$96.25

### DP Estética monochromatic tooth coloured PMMA discs

Ideal for short-term temporaries or mock-ups where strength and shade precision are valued.

**98.5mm | Shades:** A1, A2, A3, B1, C1, D2, BL1  
16mm: \$61.25 | 20mm: \$63.50

### DP Estética pink veined PMMA discs

For milling complete denture bases.

**98.5mm | Shade:** Pink veined, Original veined  
25mm: \$71.65 | 30mm: \$75.25

### DP Estética clear PMMA discs

Clear like glass. For milling hard splints and nightguards as well as orthodontic appliances. Can also be used as part of casting or pressing processes as the material burns out cleanly.

**98.5mm | Shade:** Clear  
16mm: \$51.75



unidesa · *odi*



**1-888-344-4424**  
[www.denplus.com](http://www.denplus.com)