

## VITA VIONIC® DENT DISC multiColor

Elevate your digital denture quality with **VITA VIONIC® DENT DISC multiColor**, a multichromatic composite for premium full and partial denture teeth. Used with the VIONIC DIGITAL Library and/or independent of hardware/software architecture. Offering flexible use, efficient processing and esthetics without compromise, it provides a perfect shade match to the VITA shade standards. Its composite material offers high abrasion resistance and enhanced polishing characteristics, ensuring durability comparable to conventional VITA denture teeth. With VITA, you can always expect high-quality, digitally fabricated dentures that enhance efficiency, flexibility and esthetics throughout your denture journey.




**Superior esthetics** with shade gradient for a more natural effect

**Abrasion resistance** provides high quality esthetics with improved polishing characteristics

**Flexibility** for safe, simple processing. Compatible with conventional milling systems

### Proven Quality

Wear volume\*



\* Source: Regensburg University Hospital, Polyclinic for Dental Prosthetics, Report No. 2302\_02, 4/2023

# Material Composition

VITA VIONIC® DENT DISC multiColor is an MRP composite disc comprised of inorganically filled PMMA material that matches the quality of VITA premium teeth.

- PMMA Pearls
- Silanated Filler Material
- Cross-Linked Matrix without MMA



## Integrated Shade Gradient

VITA VIONIC® DENT DISC multiColor replicates natural dentition with a shade gradient covering four levels of chromatic intensity. The impressive esthetics are immediately visible, thanks to the effective polishing characteristics.

### Available in the following shades:

Part Number	Product Description
ADA19820	VITA VIONIC® DENT DISC multiColor, A1, Ø 98.4 x h 20 mm, 1 pc.
ADA29820	VITA VIONIC® DENT DISC multiColor, A2, Ø 98.4 x h 20 mm, 1 pc.
ADA39820	VITA VIONIC® DENT DISC multiColor, A3, Ø 98.4 x h 20 mm, 1 pc.
ADA359820	VITA VIONIC® DENT DISC multiColor, A3.5, Ø 98.4 x h 20 mm, 1 pc.
ADA49820	VITA VIONIC® DENT DISC multiColor, A4, Ø 98.4 x h 20 mm, 1 pc.
ADB19820	VITA VIONIC® DENT DISC multiColor, B1, Ø 98.4 x h 20 mm, 1 pc.
ADB29820	VITA VIONIC® DENT DISC multiColor, B2, Ø 98.4 x h 20 mm, 1 pc.
ADB39820	VITA VIONIC® DENT DISC multiColor, B3, Ø 98.4 x h 20 mm, 1 pc.
ADB49820	VITA VIONIC® DENT DISC multiColor, B4, Ø 98.4 x h 20 mm, 1 pc.
ADC19820	VITA VIONIC® DENT DISC multiColor, C1, Ø 98.4 x h 20 mm, 1 pc.
ADC29820	VITA VIONIC® DENT DISC multiColor, C2, Ø 98.4 x h 20 mm, 1 pc.
ADC39820	VITA VIONIC® DENT DISC multiColor, C3, Ø 98.4 x h 20 mm, 1 pc.
ADC49820	VITA VIONIC® DENT DISC multiColor, C4, Ø 98.4 x h 20 mm, 1 pc.
ADD29820	VITA VIONIC® DENT DISC multiColor, D2, Ø 98.4 x h 20 mm, 1 pc.
ADD39820	VITA VIONIC® DENT DISC multiColor, D3, Ø 98.4 x h 20 mm, 1 pc.
ADD49820	VITA VIONIC® DENT DISC multiColor, D4, Ø 98.4 x h 20 mm, 1 pc.
AD0M19820	VITA VIONIC® DENT DISC multiColor, 0M1, Ø 98.4 x h 20 mm, 1 pc.
AD0M39820	VITA VIONIC® DENT DISC multiColor, 0M3, Ø 98.4 x h 20 mm, 1 pc.

## Fabrication of Denture Bases

### VITA VIONIC® BASE DISC HI

Shades: classic pink, classic pink translucent, soft pink, dark pink (Dark Pink available in Canada only)  
Geometries: Ø 98.5 mm / H 30 mm

