






## Upper Anteriors

 <p>ovoid</p>	<p><b>040</b></p> <p>40.0</p>  <p>10.0</p> <p>7.5</p>	<p><b>044</b></p> <p>43.9</p>  <p>10.2</p> <p>8.1</p>	<p><b>049</b></p> <p>49.2</p>  <p>11.8</p> <p>9.0</p>	
 <p>triangular</p>	<p><b>T41</b></p> <p>40.5</p>  <p>9.4</p> <p>7.5</p>	<p><b>T43</b></p> <p>43.3</p>  <p>10.3</p> <p>8.2</p>	<p><b>T46</b></p> <p>46.1</p>  <p>10.2</p> <p>8.5</p>	<p><b>T49</b></p> <p>48.8</p>  <p>11.8</p> <p>9.0</p>
 <p>rectangular</p>	<p><b>R42</b></p> <p>42.1</p>  <p>9.3</p> <p>7.9</p>	<p><b>R45</b></p> <p>44.5</p>  <p>10.3</p> <p>8.2</p>	<p><b>R47</b></p> <p>46.6</p>  <p>11.0</p> <p>8.5</p>	
 <p>square</p>	<p><b>S47</b></p> <p>46.8</p>  <p>10.8</p> <p>8.9</p>	<p><b>S50</b></p> <p>49.7</p>  <p>12.0</p> <p>9.2</p>		

\* Refer to the Available Shades chart on page 3 for the VITA SYSTEM 3D-MASTER equivalent of VITA classical shades.

## Lower Anteriors



VITA MFT Anterior	UJ / No.  mm	<b>O40</b> 40.0	<b>O44</b> 43.9	<b>O49</b> 49.2	<b>T41</b> 40.5	<b>T43</b> 43.3	<b>T46</b> 46.1	<b>T49</b> 48.8	<b>R42</b> 42.1	<b>R45</b> 44.5	<b>R47</b> 46.6	<b>S47</b> 46.8	<b>S50</b> 49.7
	LJ / No.  mm	<b>L33</b> 32.6	<b>L34</b> 34.0	<b>L40</b> 39.5	<b>L33</b> 32.6	<b>L34</b> 34.0	<b>L37</b> 36.8	<b>L40</b> 39.5	<b>L34</b> 34.0	<b>L37</b> 36.8	<b>L37</b> 36.8	<b>L40</b> 39.5	<b>L40</b> 39.5
VITA MFT Posterior	UJ / No.  mm	<b>PU29</b> 28.8	<b>PU29</b> 28.8	<b>PU33</b> 33.0	<b>PU29</b> 28.8	<b>PU31</b> 30.6	<b>PU31</b> 30.6	<b>PU33</b> 33.0	<b>PU29</b> 28.8	<b>PU31</b> 30.6	<b>PU33</b> 33.0	<b>PU33</b> 33.0	<b>PU33</b> 33.0
	LJ / No.  mm	<b>PL29</b> 29.6	<b>PL29</b> 29.6	<b>PL33</b> 33.6	<b>PL29</b> 29.6	<b>PL31</b> 31.1	<b>PL31</b> 31.1	<b>PL33</b> 33.6	<b>PL29</b> 29.6	<b>PL31</b> 31.1	<b>PL33</b> 33.6	<b>PL33</b> 33.6	<b>PL33</b> 33.6

The information contained in the combination table is only intended as a guide, and can vary depending on the patient's individual situation. \*long (O47 or L37 with longer tooth neck)

VITAPAN  
EXCELL

VITAPAN  
PLUS

VITAPAN  
LINGOFORM

VITA  
PHYSIODENS

VITAPAN

VITAPAN  
CUSPIFORM

VITAPAN  
SYNOFORM

VITA  
MFT

VITA VIONIC  
VIGO

VITA VIONIC  
FRAME

## Posterior



Available shades	
VITA classical A1-D4® *	VITA SYSTEM 3D-MASTER®
	0M1
	0M3
A1	2M1
~A2	2M2
A3	3M2
~A3,5	3R2.5
~A4	4M2
~B1	1M1
~B2	2L1.5
~B3	3L2.5
~B4	3M3
C1	–
~C2	3L1.5
C3	4L1.5
~C4	5M1
~D2	3M1
D3	–
~D4	3L1.5

\* The VITA MFT classical shades specified with "~" are not available at this time. Use the closest VITA SYSTEM 3D-MASTER equivalent shades (right column of the Available shades chart).