


QUICK MASTER FACEBOW



- 1 — **AF 201D** Right Arm
- 2 — **AF 201G** Left Arm
- 3 — **AF 202** T-bar
- 4 — **AF 203** Articulation Axis
- 5 — **AF 204** Articulation Axis Washer
- 6 — **AF 205** Arm Locking Stud
- 7 — **AF 206** Arm Locking Button

- 8 — **AF 207** Disposable Ear Piece
- 9 — **AF 208** Assembly Rod to Bow Connector
- 10 — **AF 208/B** Fork Assembly Foot
- 11 — **AF 209** Connector Locking Button
- 12 — **AF 210** Fork Support Shaft
- 13 — **AF 211** Assembly Rod Locking Connector

- 14 — **AF 212** Clamping Washer
- 15 — **AF 213** Locking Button for AF 211
- 16 — **AF 214** Fork Clamp
- 17 — **AF 215** Clamping Sleeve
- 18 — **AF 216** Inner Locking Screw
- 19 — **AF 217** Stainless Steel Bite Fork
- 20 — **AF 220/A** Handle for Disposable Bite Fork

- 21 — **AF 220** Disposable Bite Fork
- 22 — **AF 222** Nasal Assembly to T-bar Connector
- 23 — **AF 223** Locking Thumb Wheel
- 24 — **AF 224** Nasal Guide
- 25 — **AF 224/A** Nasal Guide Support Rod
- 26 — **AF 225** O-ring
- AF 800** Fork Assembly