


Lab 25 Articulator



- 1 — **QM 1001M** Lower Element
- 2 — **QM 1002** Condyle Ball
- 3 — **QM 1003** Condyle Ball Joint Set Screw
- 4 — **QM 900/925** Upper Element Lab 25
- 5 — **QM 1006** Incisal Pin Holder

- 6 — **QM 1007** Incisal Pin Locking Screw
- 7 — **QM 1008** Incisal Pin
- 8 — **QM 1009** Horizontal Incisal Table
- 9 — **QM 1010** Incisal Table and Incisal Pin Holder Racking Nut

- 10 — **QM 1011M** Upper Element Support Rod (table support)
- 13 — **QM 1400** Mounting Plate
- 14 — **QM 1402** Mounting Plate Locking Push Pin
- 16 — **QM 1102** Protrusive Lock Screw

- 17 — **QM 1103** Condyle Screw
- 18 — **QM 1104** Condyle Spring
- 19 — **QM 1105L** Condyle Saddle