


KENSON F35

F35

Introducing F35, the new extra large mould size in Kenson F-10 Posterior Teeth. This natural looking, shallow cusps posterior is wider and longer than the popular F33 mould size, making F35 ideal for partial dentures. Available in 16 matching VITA shades and whiter white shade 51.



KENSON XL & X36 TEETH

XL Teeth

Upper central incisors only...over two millimeters wider and longer*. Available in all Kenson shades.

The new Kenson XL mould saves you the trouble and expense of making teeth from scratch or reaching for more expensive brands when you have a large anterior space to fill in an upper complete or partial denture



X36 Teeth


Upper bicuspids only... Kenson mould X36 is an ideal space filler for partial dentures; can also be used as a molar.



KENSON RESIN VITA SHADES


Kenson Teeth, the professional's choice for 50 years, are now available in 16 matching VITA shades that are True to Life.

Myerson continues to offer Kenson Teeth in the 10 original Kenson shades: 55, 61, 62, 65, 66, 67, 69, 77, 81, 87, plus New whiter white shade 51.



myerson Part Number T1SGV

KENSON POSTERIOR KEY



Kenson Resin COMBINATION CHART

| MAXILLARY ANTERIOR MOULD | ARTICULATES WITH MANDIBULAR MOULD | COMBINES WITH 0° POSTERIOR MOULD | COMBINES WITH F-10 POSTERIOR MOULD | COMBINES WITH 20° POSTERIOR MOULD | COMBINES WITH 30° POSTERIOR MOULD |
|--------------------------|-----------------------------------|----------------------------------|------------------------------------|-----------------------------------|-----------------------------------|
| 14 | 43 | 29 | F29 | 29M or 29L | S30•L30•31 |
| 15 | 47 | 31 | F31 | 31M or 31L | S32•L32•33 |
| 16 | 49 | 31 | F31 | 31M or 31L | S32•L32•33 |
| 17 | 49 | 33 | F33 | 33M | 33•X34 |
| 27 | 39 | 29 | F29 | 29M or 29L | S30•L30•31 |
| 28 | 46 | 31 | F31 | 31M or 31L | S32•L32•33 |
| 29 | 48 | 31 | F31 | 31M or 31L | S32•L32•33 |
| P5 | 46 | 29 | F29 | 29M or 29L | S30•L30•31 |
| P6 | 48 | 31 | F31 | 31M or 31L | S32•L32•33 |
| P7 | 49 | 33 | F33 | 33M | 33•X34 |
| S5 | 53 | 33 | F33 | 33M | 33•X34 |
| S6 | 53 | 33 | F33 | 33M | 33•X34 |
| T3 | 43 | 29 | F29 | 29M or 29L | S30•L30•31 |
| T4 | 46 | 29 | F29 | 29M or 29L | S30•L30•31 |
| T33 | 47 | 29 | F29 | 29M or 29L | S30•L30•31 |
| T34 | 47 | 31 | F31 | 31M or 31L | S32•L32•33 |
| T5 | 48 | 31 | F31 | 31M or 31L | S32•L32•33 |
| T6 | 48 | 33 | F33 | 33M | L32•33 |

*0 and "BLOCK" Posterior teeth are interchangeable



MANUFACTURER
myerson company limited
3 trinity avenue
laventille, trinidad & tobago
phone 868 623 1007

DISTRIBUTED BY
myerson LLC
5106 north ravenwood
chicago il 60640-2713
usa
phone 800 423 2683

EUROPEAN REPRESENTATIVE
myerson limited
34 railway approach
station road
harrow middx HA3 5AA
uk
phone +44 0 20 8863 9044

Make great removables with myerson teeth & state-of-the-art myerson thermoplastic materials

myerson DuraFlex flexible partial dentures



myerson VisiClear clear clasps & frameworks



myerson DurAcetal durable cosmetic appliances



Start getting a better night's sleep

Obstructive sleep apnea & snoring can be treated safely & effectively with the ema® oral appliance from myerson



Kenson® Teeth

The Professional Choice Since 1955

KENSON RESIN ANTERIOR TEETH • MAXILLARY


| | | | | | |
|---------|--|-----------------------------|----------|--|-----------------------------|
| 14 △ | | 43.1 9.9 7.7 (43) | P7 □ | | 48.2 10.2 8.8 (49) |
| 15 □ | | 45.0 10.6 8.1 (47) | S5 ○ | | 51.5 12.1 9.4 (53) |
| 16 ○ | | 46.8 11.0 8.8 (49) | S6 □ | | 51.4 10.5 9.5 (53) |
| 17 ○ | | 49.0 11.3 8.8 (49) | T3 △ | | 41.7 9.7 7.5 (43) |
| 27 △ | | 41.5 8.7 7.4 (39) | T4 △ | | 44.6 9.8 8.4 (46) |
| 28 □ | | 45.4 9.6 8.2 (46) | T33 △ | | 44.0 10.2 8.3 (47) |
| 29 □ | | 46.9 10.3 8.5 (48) | T34 △ | | 45.6 10.8 8.3 (47) |
| P5 ○ | | 43.1 9.4 8.0 (46) | T5 □ | | 47.0 9.2 9.1 (48) |
| P6 □ | | 46.2 10.1 8.4 (48) | T6 □ | | 47.9 10.2 8.5 (48) |

TOOTH TYPE



KENSON RESIN ANTERIOR TEETH • MANDIBULAR

| | | | | | |
|----|--|--------------------|---------|--|----------------------------|
| 39 | | 31.0 8.5 4.4 | 49 | | 36.4 9.0 5.5 |
| 43 | | 32.6 9.1 4.9 | 50 | | 37.0 11.3 5.6 |
| 46 | | 33.7 8.5 5.0 | 52 | | 37.0 12.3 5.5 |
| 47 | | 33.6 9.9 4.9 | 53 | | 40.0 10.5 5.9 |
| 48 | | 35.6 8.6 5.4 | 14 △ | | 43.1 9.9 7.7 (43) |



KENSON RESIN • 0° POSTERIOR TEETH

| | | |
|----|--|------------------------|
| 29 | | 29.0 30.5 7.10 7.85 |
| 31 | | 30.5 32.3 7.55 8.50 |
| 33 | | 32.8 34.9 8.02 8.74 |

KENSON RESIN • F10 POSTERIOR TEETH

| | | |
|-----|--|------------------------|
| F29 | | 28.6 31.6 8.36 7.48 |
| F31 | | 31.0 31.8 9.19 8.37 |
| F33 | | 32.8 34.3 9.62 9.09 |
| F35 | | 36.5 37.5 10.8 9.62 |

KENSON RESIN • 30° POSTERIOR TEETH

| | | | | | |
|-----|--|-------------------------|-----|--|------------------------|
| S30 | | 28.5 29.4 9.12 8.53 | L32 | | 31.2 31.1 9.15 9.15 |
| L30 | | 28.8 30.2 9.10 9.15 | 33 | | 31.8 35.7 9.60 9.82 |
| 31 | | 31.3 33.9 9.50 10.15 | X34 | | 33.0 35.8 9.94 9.74 |
| S32 | | 31.6 32.0 9.54 9.09 | X36 | | 15.7 *Width of 2 |

KENSON RESIN • BLOCK POSTERIOR TEETH

| | | |
|----|--|------------------------|
| 29 | | 28.7 30.4 7.10 7.85 |
| 31 | | 30.5 32.4 7.30 8.50 |
| 33 | | 33.0 34.5 8.00 8.70 |

KENSON RESIN • 20° POSTERIOR TEETH

| | | |
|-------|--|--|
| 29M/L | | 29.0-31.6/29.0-31.4 8.10-8.85/8.85-8.70 |
| 31M/L | | 31.4-31.6/31.6-34.8 8.10-9.10/9.55-9.25 |
| 33M | | 33.1 35.8 8.35 9.40 |