

D I G I T A L



DIGITAL SOLUTIONS. SUPPORTED.

DenPlus proudly offers a range of solutions for scanning, 3D printing and milling from some of the global leaders in digital dentistry. We've invested in and committed to providing technical support and training for the solutions we sell, ensuring that when dental professionals invest in digital products they implement the technologies quickly, successfully and efficiently.

DETAX



keystone[®]
i n d u s t r i e s

SHINING 3D
DENTAL



VITA[®]

1-888-344-4424
www.denplus.com

myerson.

PRINT RESIN

MODEL

DenPlus 3D DLP	– beige, ivory, grey, light grey, white	1kg
Keystone KeyModel Ultra	– ivory, grey	1kg
Phrozen Aqua 4K	– ivory, grey	1kg
Phrozen Water Washable	– ivory	1kg



SPLINT

Detax Freeprint Splintmaster Flex
 Detax Freeprint Splintmaster Taff
 Keystone KeySplint Soft
 Keystone KeySplint Hard

– clear	1kg
– clear	1kg
– tinted and clear	1kg
– tinted and clear	1kg

TOOTH

Detax Freeprint Crown	– A1, A2, A3	500g and 1kg
Detax Freeprint Temp	– A1, A2, A3	500g and 1kg
Rodin Titan	– A1, A2, A3, B1, C1, 0M1, 0M3	600g
Rodin Sculpture 2.0	– A1, A2, A3, B1, C1, 0M1	600g



SURGICAL GUIDE

Keystone KeyGuide	– tinted	1kg
-------------------	----------	-----

DENTURE - FLEXIBLE

Detax Freeprint Denture Flex
 – pink, clear
 – 1kg

DENTURE

Detax Freeprint Denture	– pink	1kg
Detax Freeprint Denture Impact	– pink	1kg



PRINT RESIN



1-888-344-4424

TRY-IN

Keystone KeyDenture Try-in – A1, B1, BL4 500g and 1 kg
Rodin Try-in 1 kg



MOUTHGUARD

Keystone KeyGuard – white and black 1 kg

GINGIVAL MASK

Keystone KeyMask – pink 500g
Phrozen Gingival Mask – pink 500g



TRAY

Detax Freeprint Tray 2.0
Keystone KeyTray

– green 1 kg
– light purple 1 kg

CASTABLE

DenPlus 3D DLP casting – green 1 kg
Keystone KeyCast – purple 1 kg

www.denplus.com

*Visit our web site for
current pricing and
promotions.*

3D PRINTERS

SHINING 3D
DENTAL

ACCUFAB L4D

- **4K resolution**
- Add printing capacity affordably and without compromise
- Small and large build plates/resin vats for maximum efficiency and flexibility
- Free slicing software with no annual license fees
- Free add-on modules for model, bite splints and crown design
- Resin vat lifespan controlled by user - not via RFID chip
- Large selection of pre-validated resin profiles as well as ability to create custom profiles

ACCUFAB CEL

- **6K resolution**
- Add printing capacity affordably and without compromise
- Small and large build plates/resin vats for maximum efficiency and flexibility
- Free slicing software with no annual license fees
- Free add-on modules for model, bite splints and crown design
- Resin vat lifespan controlled by user - not via RFID chip
- Large selection of pre-validated resin profiles as well as ability to create custom profiles
- Print chamber temperature control/heater



Available Q4 2025



ACCUFAB F1

- DLP projection system
- 1920 x 1080 resolution
- Accuracy +/- 25um
- Dual heating technology significantly reduces printing time and improves printing success rate
- Intelligent sensing: the printing platform, resin tank, and optical machine feature automatic detection to monitor missing components or residues in real-time
- Free slicing software with no annual license fees
- Large selection of pre-validated resin profiles as well as ability to create custom profiles

www.denplus.com

Visit our web site for current pricing and promotions.

3D PRINTING ACCESSORIES



1-888-344-4424

LIGHT CURING

SHINING 3D FABCURE 2

- LED Power 70W
- Large curing chamber
- Built-in heater
- Pre-set curing times by resin type/brand
- Programmable for custom strategies



SHINING 3D FABCURE 1S

- LED Power 30W
- compact footprint where space is limited

PHROZEN CURE+

- Flash curing
- Heated enclosure
- Broad wavelength range to ensure optimal curing for a range of materials
- Programmable for custom strategies
- LED power 150W



WASHING

PHROZEN WASH +

- Two compartments for simplified two-stage cleaning of printed parts
- Programmable washing cycles
- Designed to accommodate build-plate from XL 4K printers or place printed parts into reservoirs



SHINING 3D FABWASH

- Pre-programmed 'wash and dry' cycles
- Single compartment with 9.5L capacity
- Designed to accommodate build-plate from Accufab L4D, CEL and D1S printers or place printed parts into reservoir
- Automatic part removal from Accufab L4D, CEL and D1S build plates

SCANNERS

SHINING 3D
DENTAL

DESKTOP SCANNER

SHINING 3D DS-EX PRO H

- **Scan accuracy 8 μ m**
- **Dual 5.0 MP cameras**
- STL, OBJ and PLY scan file formats
- Free scanning software with no annual license fees
- Seamless integration with exocad CAD
- Multiple scan modes
- Supports full-size articulators and impression scanning



www.denplus.com

*Visit our web site for
current pricing and
promotions.*



SHINING 3D DS-EX PRO

- **Scan accuracy 10 μ m**
- **Dual 1.3 MP cameras**
- STL, OBJ and PLY scan file formats
- Free scanning software with no annual license fees
- Seamless integration with exocad CAD
- Multiple scan modes
- Supports full-size articulators and impression scanning

SCANNERS



1-888-344-4424

INTRA-ORAL SCANNER

SHINING 3D AORALSCAN 3

- **Wired OR Wireless**
- STL, OBJ and PLY scan file formats
- Free scanning software with no annual license fees
- Two scanner tips sizes. Tips autoclavable up to 100 times.
- Software includes free model creation and temp design modules



SHINING 3D AORALSCAN ELITE IOS

- **Wired OR Wireless**

The Aoralscan Elite is the world's first device to integrate both intraoral scanning and photogrammetry into a single unit. The device functions as a standard intraoral scanner (IOS) for capturing images of dentulous cases, allowing for detailed scans of teeth and oral structures. Additionally, it serves as a photogrammetry system, which is used to accurately record the position of dental implants in edentulous cases.



FACE SCANNER

SHINING 3D METISMILE

- 10 seconds to acquire complete facial scan
- Free scanning software with no annual license fees
- STL, OBJ and PLY scan file formats
- Aligned scan files easily imported into all CAD design software
- Optimized for exocad CAD



PMMA DENTURE BASE

DenPlus Estética Base (98mm) – pink veined
25mm, 30mm original veined

Keystone KeyMill (98mm) – original and
25mm, 30mm light reddish pink

Yamahachi Pink Veined (98mm) – pink veined
25mm, 30mm, 35mm



PMMA TOOTH

DenPlus Estética (98mm)
Mono: 16mm, 20mm
Shades: A1, A2, A3, B1, C1, D2, BL1

8-layer translucent:
16mm / Shades: A1, A2, A3, A3.5, A4, B1, C1, C2, C3, C4, D2, D3, D4, BL1

20mm / Shades: A1, A2, A3, B1, C2, BL1

10-layer: 20mm
Shades: A1, A2, A3, B1, C2, BL1

Yamahachi (98mm)
Mono: All 16 classic Vita shades
12mm, 14mm, 16mm, 18mm, 20mm, 25mm

Vita Vionic Dent Disc Multicolour (98mm - 20mm)
Shades: All 16 classic Vita shades, 0M1, 0M3



PMMA CLEAR

DenPlus Estética Clear (98mm):
14mm, 16mm, 18mm, 20mm, 25mm

Yamahachi Clear (98mm):
12mm, 14mm, 16mm, 18mm, 20mm, 25mm, 30mm



POLYMER

Myerson Duraflex CAD (98mm)	– medium pink, dark pink	15mm, 20mm, 25mm
Myerson Duracetal CAD (98mm)	– A1, A2, A3, A3.5, B1, C3, bleach, pink	15mm, 20mm, 25mm
Myerson Visiclear CAD (98mm)	– clear	15mm, 20mm, 25mm
Myerson Dentivera (98mm)	– tan and white	18mm, 24mm