

PHOTO-TRAY AND PHOTO-BASE

- 50 sheets per box

\$69 ~~\$86⁹⁹~~
Save
20%

3051101

3051102



MOLLOPLAST-B

detax

\$329⁶⁶ ~~\$359⁶⁶~~ 170 g
DET03001
Save
\$30



ORTHODONTIC BOXES

- 25 per box

\$19 ~~\$28⁷⁵~~
Save
33%



DENPLUS BASEPLATE WAX

\$59 ~~\$66²⁵~~ 5 lb
Save
10%



3001101

3001102

VITA VIONIC DENT MULTILAYER PMMA MILLING DISCS

VITA®

Buy **10 VITA VIONIC DENT DISC**
and get either **VITA VIONIC BOND MULTIDOSE** or **VITA VMLC PRIMER**



ESTÉTICA PMMA MILLING DISCS

- Pink and Original
- 25 mm and 30 mm

\$52 ~~\$59⁹⁹~~

Save
13%



DP Estética  unidesa · odi

- 16 mm
- 8 layers translucent – Assorted Vita Classic shades

\$82⁹⁹

Buy
5
GET ONE FREE

STERILIZERS

Now available at Denplus

ENBIO S AND PRO

- ◆ Swiss made Class B air-removal (pre-vacuum) autoclave
- ◆ Ultra fast sterilization cycle
- ◆ Dimensions S - 56 cm x 25 cm x 16 cm
- ◆ Dimensions PRO - 56.1 x 26.9 x 20.3 cm
- ◆ Small footprint – ideal when space is limited



CLAVE 8+

- ◆ Class S Steam Sterilizer with air assisted Drying Technology
- ◆ Ultra-fast Closed door drying with air assisted drying
- ◆ 5 Programmed Sterilizations Cycles: Solid/Unwrapped, Wrapped, Liquid, Textile, Prions
- ◆ 8 L chamber volume (17 cm x 32 cm)
- ◆ Capacity for three trays

\$3,995 ~~\$4,620~~



CLAVE 16+

- ◆ Class S Steam Sterilizer
 - ◆ 5 Programmed Sterilizations Cycles: Solid/Unwrapped, Wrapped, Liquid, Textile, Prions
 - ◆ Air Assisted Drying Technology for Faster Drying Times – 30 Minutes Drying
- ◆ 16 L chamber volume (22.9 cm x 36 cm)
- ◆ Capacity for three trays

\$5,899 ~~\$6,773~~



CLAVE B

- ◆ Class B steam sterilizer
- ◆ Fast Sterilization and Drying Time
- ◆ Pre/Post Vacuum System for Dynamic Air Removal
- ◆ Better Penetration and faster more effective drying
- ◆ 23 L chamber volume (24 cm x 45 cm)
- ◆ Capacity for four or five trays

\$8,999 ~~\$10,463~~

